

revvity	SonoEQ Version: 2.0.2		Revision: 22
	Title: SonoEQ User Guide		
	Proprietary & Confidential	Effective Date: 11/22/2023	Page 1 of 108

SonoEQ® User Guide

This document describes how to use the SonoEQ® software package to control Vega® ultrasound imaging systems. It includes descriptions of how to perform wide-field imaging scans, quantify features of interest from the 3D image data (e.g. a tumor), and export the results.

For Technical Support or other inquiries please visit:
<http://www.revvity.com/contact-us>

Revvity, Inc.
940 Winter Street
Waltham, Massachusetts 02451 USA
www.revvity.com

Table of Contents

1	Welcome	4
2	Background	5
2.1	Acoustic Imaging Modes	5
2.2	Ultrasound Transducers	6
2.3	Scan Dimensionality	7
3	Getting Started	9
3.1	Starting the SonoEQ Software	9
3.2	Initializing the Hardware	11
3.3	Checking the Temperature of the Heated Imaging Bays	12
3.4	Setting Study Information	13
4	Acquiring Image Data	15
4.1	Acquiring Acoustic Image Data	15
4.2	Previewing Acquired Data	24
4.3	Acquisition Settings/Presets	25
5	Analyzing Image Data	34
5.1	Loading Image Data	34
5.2	Viewing/Navigating Image Data	40
5.3	Overlaying Multiple Datasets	41
5.4	Creating Segmentations and Quantifying Data	45
5.5	Additional Utilities	76
5.6	Display Options	79
5.7	Viewing Study Info	86
5.8	Analysis Settings	87
6	Exporting Data	90
6.1	Export Statistics	90
6.2	Exporting Image Data	93
6.3	Study Report Generator	94
7	Tools	97
7.1	Data Processing	97
7.2	Volume Compressor	100
8	Instrument	101
8.1	Servicing Dashboard	101
8.2	Heat Imaging Bays	101

revvity	SonoEQ Version: 2.0.2		Revision: 22
	Title: SonoEQ User Guide		
	Proprietary & Confidential	Effective Date: 11/22/2023	Page 3 of 108

8.3	Initialize Motion Stage	101
8.4	Firmware Updater	101
9	Data Storage	104
10	Mouse and Keyboard Shortcuts	105
11	Troubleshooting	107

1 Welcome

Vega imaging platforms are controlled by SonoEQ, the **E**xploration and **Q**uantification software package. SonoEQ software comes preinstalled on (1) desktop computer that is provided with the purchase of your system, enabling both data acquisition and data analysis. This system is the only approved computer to control the hardware and perform data acquisition.



Additional installation of SonoEQ on different computer workstations for data analysis is possible depending on the purchased software license. Please contact your Vega sales representative for more information.

System Requirements:

Computers running SonoEQ need enough processing, memory, and graphics capabilities to manipulate and render 3D image data, with individual file sizes equal to or exceeding hundreds of MB.

- OS: Microsoft Windows 10 or higher (64-bit)
- Processor: Minimum of 4 CPU logical cores.
- Graphics: Dedicated GPU with support for OpenGL 3.2 or later (1 GB is recommended).
- Display: Two supported resolutions of **1920 x 1080** or **1920 x 1200**.
- Memory: Minimum RAM of **4 GB** (8 GB or more is recommended).
- Mouse: A three-button mouse with scroll wheel is recommended.
- Other: A PDF viewer (Adobe Acrobat Reader is recommended)


Legend

NOTE		A NOTE indicates important information that helps you make better use of your system.
CAUTION		A CAUTION indicates either potential damage to hardware or loss of data and tells you how to avoid the problem.

2 Background

2.1 Acoustic Imaging Modes

The instrument includes several acoustic imaging modes that vary in their function and purpose. Modes are selected by clicking on the respective **Mode** button located in the **Acquisition Tab**. The following sections give a brief background of the mode functionality and potential applications.



2.1.1 Pulse/Echo Ultrasound (B-Mode)

Pulse/echo ultrasound imaging is performed by sending high frequency ultrasonic pulses (ranging from 1 - 50 MHz) into a subject and listening for returned echoes using the same transducer. As ultrasound waves travel through the body, they encounter interfaces between tissues and get partially reflected toward the probe. The proportion of the wave that is reflected versus transmitted across the interface is related to the tissues' underlying material properties (e.g. density, speed of sound, etc.). Images are made by recording the ultrasonic echoes and assigning grayscale values to their intensities, which is referred to as brightness (B)-mode.

B-Mode imaging is highly effective at locating anatomical structures and can have very fast acquisition frame rates (up to or beyond 100 Hz). B-Mode is used for:

- Identification of anatomical landmarks and organs
- Volumetric measurements (e.g. tumor volume, kidney volume, etc.)
- Generating a background orientation image for other advanced imaging modes
- Real-time scout scanning for a feature of interest
- Measuring echogenicity of an organ (e.g. fatty liver)

2.1.2 Acoustic Angiography

Acoustic angiography (AA) is a patented contrast-enhanced ultrasound (CEUS) imaging mode, which is made possible by custom dual-frequency ultrasound transducers. In AA mode, microbubble contrast agents are excited by a low-frequency ultrasound transducer element, like those used in clinical imaging (2-4 MHz), yet acoustic echoes from microbubble contrast agents are received with a second transducer at a much higher frequency (25-30 MHz). The excited microbubbles produce a broadband "super-harmonic" response, which can be detected with the high-frequency receiver, while tissue produces negligible echoes at these frequencies and is cancelled from the image. This technique yields images of blood vessel morphology entirely unlike standard B-Mode ultrasound images and resembling X-ray angiograms, thus leading to the name "acoustic angiography".

Acoustic angiography is used for evaluating microvasculature. Additionally, it can be used for:

- Detecting areas of hyper/hypo-vascularization (e.g. tumors, avascular cysts)
- Measuring blood vessel density
- Visualizing blood vessel architecture

2.1.3 Shear Wave Elastography

Shear wave elastography (SWE) is a technique used for assessing tissue stiffness by measuring how fast a mechanical wave travels through tissue. A long-duration and high-intensity ultrasound pulse is delivered to “push” the tissue using acoustic radiation force and create a shear wave, which is tracked using successive imaging pulses. The speed of the wave can be measured using the time at which the generated wavefront reaches different lateral locations. Shear wave velocity is proportional to tissue stiffness, and if the tissue type satisfies several simplifying assumptions (isotropic, linear, elastic, incompressible), the speed can be converted into a quantitative modulus value.

The following equations are utilized by SonoEQ to convert shear wave velocity (c_s) into either shear (G) or Young’s (E) modulus, assuming a tissue density (ρ) equal to 1 g/cm³ and a Poisson’s ratio (ν) equal to 0.5:

$$G = \rho c_s^2$$

$$E = 2G(1 + \nu)$$


$$E = 3G = 3\rho c_s^2$$

Shear wave elastography can be used for:

- Assessing stiffness in tissue

2.2 Ultrasound Transducers

Vega systems can be configured with several transducers preinstalled into the robotic carriage. There are two main classes of transducers, “linear arrays” and “wobblers”, which can be selected by accessing the **Settings** and navigating to [Imaging Mode] → Scan Protocol.



2.2.1 Linear Array

A linear array transducer consists of hundreds of ultrasonic crystals arranged in a row. Groupings of elements are fired simultaneously to produce individual scan lines that can then be rapidly put together into a single 2D frame. Linear array transducers typically produce rectangular fields of view with uniform beam density throughout. Advantages of the linear array include: adjustable focal depth, multiple focal depths, angular compounding for speckle reduction, deeper imaging, and faster frame rates.

revvity	SonoEQ Version: 2.0.2		Revision: 22
	Title: SonoEQ User Guide		
	Proprietary & Confidential	Effective Date: 11/22/2023	Page 7 of 108

Linear arrays are useful for imaging:

- Orthotopic and large subcutaneous tumor sizing
- Cardiovascular dynamics
- Deep abdominal imaging (e.g. liver)

2.2.2 Wobbler


A wobbler transducer consists of a single element piston that is mechanically “wobbled” to generate a 2D image. The piston is typically geometrically curved to a single focal point. Wobbler transducers produce a trapezoidal (fan)-shaped field of view due to the divergence of the ultrasound beam with increasing depth and wobble angle. This allows for a wider field of view, but with decreased line density at depth. Wobbler imaging typically has a reduced frame rate compared to linear array imaging; however, it can create higher quality images at the focus point due to higher frequencies and isotropic aperture geometry. Additionally, wobbler transducers exclusively support acoustic angiography (AA) mode.

Wobbler transducers are useful for imaging:

- Small subcutaneous tumors
- Shallow abdominal organs
- Epidermal imaging


2.3 Scan Dimensionality

Data dimensionality is an important concept to keep in mind while collecting and analyzing data in SonoEQ. In general, imaging data can be classified by up to three spatial dimensions (axial × lateral × elevational) and one temporal dimension (time). For the remainder of the manual, we will utilize the terms “sequence” and “static” to differentiate between data captured with and without a temporal dimension, respectively. Specifically, a “static” dataset is one where there are only spatial dimensions, for example a single 2D B-Mode frame or a single 3D acoustic angiography volume. A “sequence” dataset, on the other hand, is one that includes a time dimension, such as a 2D B-Mode movie of a heart beating, where the transducer is held in one place and 2D images are captured with a real-time frame rate. Be aware that not all imaging modes or transducers support all dimensionality in both “static” and “sequence” capture types. Scan dimensionality is chosen by clicking the appropriate section in the **Acquisition Tab**.




2.3.1 M-Mode (1D) Scan

The motion (M)-Mode tab is primarily used for capturing 1D “sequence” datasets. An M-Mode is typically represented as a 2D image, where the vertical axis is the axial spatial dimension, and the horizontal axis is the time dimension.



2.3.2 2D-scan


The 2D scan tab is used for capturing “static” 2D images and 2D “sequences”. In these images, the vertical axis is the axial dimension, and the horizontal axis is the lateral (azimuth) axis. This scan type is the typical view offered by conventional ultrasound scanners.



2.3.3 3D-scan

The 3D scan tab is used for capturing “static” 3D volumes and 3D “sequences” (also referred to as 4D-Mode). 3D datasets are reconstructed by stitching together many 2D frames into a single volume. In the case of a 3D-sequence, 2D sequences are captured in a stepwise fashion, and you can optionally apply **Retrospective Phase Gating** to align the phases of the 2D movies into one cohesive 3D volume sequence. More detail about phase gating can be found in:

D. Chittajallu et al., “Image-based methods for phase estimation, gating and temporal super-resolution of cardiac ultrasound,” IEEE Trans. Biomed. Eng., vol. 66, no. 1, pp. 72–79, 2018.




revvity	SonoEQ Version: 2.0.2		Revision: 22
	Title: SonoEQ User Guide		
	Proprietary & Confidential	Effective Date: 11/22/2023	Page 9 of 108

3 Getting Started

3.1 Starting the SonoEQ Software



With the system powered on, and an animal properly sedated and placed on the imaging pad, you are ready to collect data.

- 1) To start the SonoEQ software, double click the  icon on the desktop.



NOTE:

SonoEQ.exe is installed by default under the parent folder “C:\Program Files\Revvity\SonoEQ 2.0.2” and a shortcut is placed on the desktop.

If the desktop icon has been removed, click the Windows Start menu , or Windows key on keyboard, type “SonoEQ”, and click the  SonoEQ icon.


- 2) You should see the SonoEQ Splash Screen appear to indicate the startup progress.



Upon loading the splash screen will disappear and the main SonoEQ window will appear.


The SonoEQ software is composed of two main graphical components:

- Control Panel
 - This is where options and settings are set, 3D scans are deployed, etc.
- Slice Views
 - This is where image data can be viewed.



From within the control panel, the two main modes of software operation can be selected. By default, SonoEQ launches into **Analysis** mode.

- **Acquisition mode:** Collect new imaging data
- **Analysis mode:** Quantify existing data (e.g. determine tumor volume)




3.2 Initializing the Hardware


When entering Acquisition mode for the first time from a powered down state, the instrument must be initialized. The initialization procedure involves connecting to each hardware subsystem (ultrasound engine, robotic motion controller, etc.) and running short calibration procedures for motorized components. Initialization moves all motor driven components through their full range of motion to establish home position and resets electronics and controllers.

For further details on instrument operation, see the ***Vega Hardware Manual*** provided with your device.

- 1) Click the **Image Acquisition** tab to enter Acquisition Mode.




- 2) A progress dialog will appear and the software will connect to each of the system components and initialize the motion stages to the home position (if necessary). Homing initialization takes approximately 1 min to complete.



NOTE: Each time the computer is powered down, the software must re-home the motion stages. It is recommended to leave the computer power ON for the duration of an imaging timepoint to avoid multiple initialization operations.

- 3) Upon completion of initialization, the software will enter B-mode by default, and you will see both the camera and the ultrasound images streaming in real time within the Slice Views window. The camera stream displays the surface of the imaging platform, and the ultrasound stream is the 2D data coming from the ultrasound transducer's current position.






NOTE:


The ultrasound stream can be paused/frozen using the **Freeze Probe** button located in the bottom-left of the Control Panel. Freezing the probe ceases all sound transmission from the currently active probe.

3.3 Checking the Temperature of the Heated Imaging Bays


Vega devices with heated imaging bays will have an indicator in the lower right corner of the “Image Acquisition” screen. This indicator informs you about the current state of your heated imaging bays. When the Vega is powered ON, the heater will start to warm the animal bays. Average heat-up time from room temperature (21°C) to body temperature (38°C) is ~20 min.

<u>Icon</u>	<u>State</u>	<u>Description</u>
Imaging Bays: 21C/38C 	NOT READY	Imaging bays are not at the desired temperature.
Imaging Bays: 38C/38C 	READY	Imaging bays are at the desired temperature.
Imaging Bays: 36C/--C 	OFF	Imaging bays are off and not heating.

Heated imaging bays can be turned ON and OFF using the **Instrument** menu.



Heat Imaging Bays ON




Heat Imaging Bays OFF

Upon closing acquisition, you will be prompted to choose whether you would like to keep heating the imaging bays. It is recommended you keep the imaging bays heated if you plan to continue imaging within a short period of time to avoid unnecessary heat-up time.


3.4 Setting Study Information

Study information is entered in the textboxes located in the left **Control Panel**. Data is saved to the Windows File System using the directory path specified in the corresponding textbox.

- 1) Enter a new "Study Name." "Study Directory" can be edited using the  button.




- 2) If desired, enter a new "Timepoint" name. By default, this is set to the current date.












NOTE:

Study names and timepoints are used to index your data acquired on specific days, or of specific mice over time. For instance, a Study Name might be "ChemotherapyResponseOverTime", and Timepoint name might be "Day5AfterTumorImplants"

- 3) Enter the ID's of up to three animals in the "Animal ID" textboxes. These might correspond to the animal's unique animal identification scheme used by your lab, such as "Cage202_Animal2".



- 4) (Optional) Click the **Animal Details**  button to add additional information about the study. The data entered into the table will be saved as a header in each image file, thereby serving as a permanent record of study parameters. Use the  button to propagate fields across animal columns (for example, if all three animals are the same *Strain*).

		Animal1	Animal2	Animal3	
Animal Model	<input checked="" type="checkbox"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Animal Strain	<input checked="" type="checkbox"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Group	<input checked="" type="checkbox"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Injection Time	<input checked="" type="checkbox"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Time After Injection	<input checked="" type="checkbox"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Scan Name	<input checked="" type="checkbox"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Comment	<input checked="" type="checkbox"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	

Clear Table


OK

4 Acquiring Image Data

4.1 Acquiring Acoustic Image Data

The following section describes the process for capturing acoustic (i.e. ultrasound) image data with SonoEQ. The most common use-case is presented first – capturing a 3D static B-mode volume – followed by subsections that describe how to acquire data for each mode and dimensionality type.
















As described above, nearly all imaging modes can be set to capture single “static” image, or multi-frame movie “sequences”. For more background on static and sequence captures, the reader is referred to section 2.3.



NOTE: For all Image Acquisitions, a thumbnail image of the scan region is automatically saved along with the scan data.

4.1.1 Acquiring a 3D B-mode/AA Volume


- 1) After an animal ID has been entered, the corresponding column in the Scan Plan Table will become active, and a rectangular region of interest (ROI) will appear in the Live Camera Stream. Clicking the eyeball button toggles the visibility of the scan planning object. If the scan planning object is hidden, that bay will not be scanned during acquisition.

M-Mode	2D Scan	3D Scan	
	Bay 1	Bay 2	Bay 3
Run Heart Finder		 Find Heart	
Visibility			
ROI W/L (mm)	30.00  40.00 	30.00  40.00 	30.00  40.00 
	Default	Default	Default
	Passes: 3 	Passes: 3 	Passes: 3 
Duration (static/sequence)	0m22s / N/A		



NOTE: The scan plan table offers useful information about the size of the scan region (W x L) as well as the estimated time required for each scan to be completed. Larger scan regions will increase the total duration of the data acquisition.

- 2) Select the region to scan in the camera view using the rectangular scan planning objects. Clicking and dragging on the green handle in the center re-positions the entire scan region, while clicking/dragging the red and blue handles adjusts the bounds.




NOTE: The location of the ultrasound probe is indicated on the camera view and can be manipulated with the motion control pad for a quick glance into the animal.





- 3) Once you are satisfied with the scan region placement, click **Acquire Static** to begin the acquisition.

Acquire Static

- 4) Once the data is captured and processing is finished, clicking the **Load**  button under the Animal ID textbox will load the scan for previewing (see section 4.2).

4.1.2 Acquiring a 2D B-mode Image



- 1) To collect 2D scans, navigate to the **2D Scan** tab in the Control Panel.
- 2) The transducer can be moved to the desired location using the directional arrows in the motion control pad. The circular Go-To-Position button in the center allows you to send the transducer to a selected location by clicking the button and then clicking on the slice viewers (webcam or ultrasound live stream, or on a loaded volume).
- 3) A 2D static image or sequence can be acquired using either the **Acquire Static** or **Acquire Sequence** button, respectively. 2D sequences can be captured to assess cardiac function using the same protocol as acquiring an M-Mode, and the “2D Scan” tab also includes the Heart Finder row to acquire heart location images.
- 4) Once the data is captured and processing is finished, clicking the **Load**  button under the Animal ID textbox will load the scan for previewing (see section 4.2).




NOTE:

Retrospective cardiac phase gating can be enabled for 2D sequences (see section 4.3.2).


4.1.3 Acquiring a 1D M-Mode Sequence

M-Mode imaging produces high frame rate sequences to characterize moving anatomy, such as left ventricular motion of the heart. The image that is displayed represents one line of the B-Mode image plotted over time.



Animal Positioning: For optimal M-Mode cardiac images, the animal should be in the prone position and be angled by approximately 30 degrees to the right of vertical to obtain a short axis view of the heart.

- 1) To acquire M-Mode data, click the tab labeled M-Mode in the Acquisition panel. The software will switch into a 3-slice view showing real-time streams of both B-Mode and M-Mode. The position of the single M-Mode line is indicated by a green vertical line shown in the B-Mode stream.






NOTE: M-Mode images can only be acquired with a linear array transducer.

- 2) (Optional) If desired, you can run the Heart Finder algorithm to automatically find the heart in 3D space. This feature is useful to consistently place the M-Mode interrogation line in the appropriate position when quantifying cardiac parameters such as Ejection Fraction (EF). To

revvity	SonoEQ Version: 2.0.2		Revision: 22
	Title: SonoEQ User Guide		
	Proprietary & Confidential	Effective Date: 11/22/2023	Page 19 of 108

enable the Heart Finder workflow, first click the **Place ROI** button located in the table. An ROI will appear (like those shown in 3D Image mode). The button text will change to **Find Heart**, and upon another click, the region of interest will be scanned, and a heart location image will appear.



- 3) (Optional) After the heart has been localized, click the center button of the motion control pad, and then click the location on the Heart Location image you wish the transducer to move to.



- 4) Once you are satisfied with the scan region placement, click **Acquire Sequence** to begin the M-Mode acquisition. Sequence duration can be set in the Settings menu (see section 4.3.1) and is set to a default of 2 sec. The corresponding B-Mode is also saved.





- 5) Once the data is captured and processing is finished, clicking the **Load** button (indicated by a blue downward arrow) under the Animal ID textbox will load the scan for previewing (see section 4.2).

4.1.4 Acquiring a 3D B-mode Sequence


The only application currently supported for 3D sequences is for cardiac imaging, and the linear array is the only supported transducer. Data from a 3D sequence will automatically be

processed as a cardiac data set and will result in a 3D sequence of a heart beating over a single cycle.

- 1) In the 3D tab, the **Place ROI** button in the Heart Finder row is automatically set to **Find Heart**. Clicking the button will run a scan to produce the heart location image using the existing ROI. Previous heart finder scans, like those collected for M-Mode acquisitions, can be loaded through the **Load**  button menu and used for positioning of acquisition ROI (discussed further in section 4.2).



- 2) Change the ROI width to a single pass (12.8 mm) and a length of about 10 mm. Move the ROI to be centered around the bottom of the heart location image, ensuring the ROI includes the apex of the heart.



- 3) Select the desired **Sequence Duration** (see section 4.3.1). The phase correction algorithm is optimized for a duration of 2 seconds.



NOTE: A longer duration may increase the processing time of phase correction.

- 4) Click **Acquire Sequence** to start the scan.




- 5) Once the data is captured and processing is finished, clicking the **Load** button under the Animal ID textbox will load the scan for previewing (see section 4.2).


4.1.5 Acquiring a 2D SWE Image


The shear wave elastography (SWE) mode consists of a 2D B-Mode stream, where a smaller region must be defined for capturing a stiffness image. In SWE mode, only 2D and 3D images can be acquired so the **M-Mode** tab will be disabled. To collect 2D SWE images, navigate to the **2D Scan** tab in the Control Panel.

- 1) The transducer can be moved with the motion control pad (orange box) to center the green analysis window (white arrow) on the tissue of interest.




- 2) The axial and lateral bounds of the analysis window can then be adjusted to further focus on a region by dragging the green or red interaction handles on the analysis window box. The green analysis box cannot extend beyond the blue box, so if the desired imaging target lies outside of the blue box the transducer needs to be moved or the animal repositioned.
- 3) Click the **Acquire Static** button to capture a 2D SWE volume as well as a reference B-Mode.




- 4) Once the data is captured and processing is finished, clicking the **Load**  button under the Animal ID textbox will load the scan for previewing (see section 4.2).

4.1.6 Acquiring a 3D SWE Volume

The workflow for acquiring a 3D SWE volume is just like acquiring a 3D B-Mode volume, with the addition of adjusting the SWE analysis window. At this time, 3D SWE acquisition can only be acquired as a Stepped Scan due to the limited frame rate of the SWE mode (for more details on Step Mode settings see section 4.3.1).



- 1) A scan region can be selecting by moving and resizing the ROI on the coronal (green) slice either on the live camera view or a loaded volume.
- 2) The axial bounds of the analysis window can then be adjusted to further focus on a region by dragging the green or red interaction handles on the analysis window box.
- 3) Click the **Acquire Static** button to capture a 3D SWE volume as well as a reference 3D B-Mode.
- 4) Data can be loaded into the scene for preview by clicking the **Load**  button under the Animal ID textbox (see section 4.2).




NOTE:

The green analysis region on the ultrasound live stream cannot be made wider than the blue box, but the lateral extent of the scan is determined by the ROI on the coronal (green) slice, and not the width of the blue box.

4.2 Previewing Acquired Data

- 1) After clicking **Acquire Static** or **Acquire Sequence** for the different imaging modes and dimensionalities, a progress bar will indicate how long until the scan is complete, and various icons will appear under the Animal ID textboxes to indicate when a volume has been reconstructed and is ready for loading/viewing.





Empty cell indicates no data has been captured for that animal yet.



Wheel icon indicates waiting actions such as when raw volume data is currently being reconstructed or acquisitions are being loaded into the scene.



Load icon indicates that volumes have been reconstructed successfully, saved to disk, and can be loaded in to preview.


- 2) Once all scanning is complete, you can either click on the **Load**  button under the Animal ID textbox to bring up the Analysis view, or the **Next**  button to clear the animals' info and begin scanning the next set of animals.



NOTE:

Clicking the **Next** button will clear *ALL* animal IDs from the textboxes and remove loaded volumes from the scene. To clear individual animals, use the drop-down arrow as shown below. Note that clicking the **Next** button does not delete data, it simply clears the textboxes.

- 3) The arrow on the side of **Load** button is a menu that contains all previously acquired datasets for that study name/timepoint/animal ID combination.



NOTE:

Default save location: %User Profile%\Documents\SonoEQ\Study_Name\Timepoint\Animal_ID\B-Mode_3D Scan_YYYY_MM_DD_HH_MM_SS\

4.3 Acquisition Settings/Presets


Imaging settings can be adjusted to improve image quality for different applications and obtain various kinds of images. Acquisition Settings are divided into **Scan Protocol** and **Post-Scan** settings for each imaging mode. Furthermore, each imaging mode has factory-defined presets that update multiple settings in the **Scan Protocol** and **Post-Scan** panels. Specific settings can still be customized, but clicking a scan preset button will return specified settings to their preset defined values.

To access the settings, either click on the **Settings Tray (>>)** button that appears next to **Control Panel**, press **Ctrl + “2”**, or click on **Edit → Settings**.

4.3.1 Scan Protocol Settings

Settings in the Scan Protocol panel control all aspects of the image acquisition and are grouped into **Transducer** settings and **Scan-Protocol** settings.

Transducer settings – This group of settings control the performance of the ultrasound transducers.



- **Imaging Depth** – Total range/extent of imaging window.
- **Time-Gain Compensation (TGC)** – Gain that can be changed to change image brightness. The overall brightness can be updated using the Overall TGC slider which will adjust the individual areas accordingly.

- **Individual TGC** – Sliders to adjust the image brightness as different depths to account for signal attenuation.
 - **Apply Default TGC** – Reverts the TGC sliders to the factory settings.
- **Focal Depth** – Depth at which ultrasound is focused and the image is the sharpest.
 - **Multi-Focal Zone** – Up to four focal zones can be activated at once for improved image quality through depth. Each additional focal zone that is added will reduce the imaging frame rate by the corresponding number of zones. E.g. if the frame rate with 1 focal zone is 100 Hz, 2 focal zones will reduce to 50 Hz, 3 focal zones to 33 Hz, etc.
- **Transmit Frequency** – Frequency of output ultrasound pulses in units of MHz. Higher frequencies have better resolution, but worse depth of penetration.
- **Frame Rate** – Rate at which 2D images are acquired.
- **Frozen entering imaging mode** – When checked, the probe will be frozen when the specified Mode is entered. E.g. if the acquisition is in B-Mode and switches to AA-Mode, if checked, the probe will automatically freeze prior to switching into AA. This setting can prevent unintentional transmission of low frequency sound waves when toggling modes.
- Wobbler only:
 - **Frame Acquisition** – Number of directions over which to acquire data on the wobble path. Frame rate is halved when using “unidirectional” mode.
 - **Jitter Compensation** – Value to compensate for any misalignment between images taken from each direction of the wobble (see “bidirectional” mode).
- Linear Array only:
 - **Spatial Compounding** – Averaging of frames to reduce speckle. Sound waves for frames are “steered” and combined to provide images from different angles.
 - **SC Count** – Number of frames to be averaged.
 - **SC Angle Delta** – difference between the steering angles for each frame.
 - **Collect RF Data** - Collect raw radiofrequency ultrasound data.

Shear Wave Elastography Transducer Settings – This group of settings controls the linear array settings for shear wave elastography acquisitions.

Transducer: Linear Array

*Selecting a different transducer will change displayed transducer settings and some scan plan settings

▼ Transducer Settings

▼ Advanced Settings

Push Offset: 0.8 mm to the Left of the tracking region

► Scan Plan Settings

- **Push Offset** – offset between the lateral location of the push pulse and tracking region.
 - The push pulse can be transmitted to the right or left of the tracking region.

Acoustic Scan-Plan Settings – This group of settings affects 3D ultrasound acquisitions and sets the shape/timing of the acquisition.

▼ Scan Plan Settings

Frame Spacing:

Scan Speed:

Default ROI Width:

Default ROI Length:

Remember ROI Size and Position: ☐

Step Mode:

Frames Per Position:

Sequence Duration:

Current Scan Plan:

▼ Advanced Timing Settings

Minimum Δt Between Stepped Positions:

Minimum Δt Between Starts of Passes:

Repetitions:

Minimum Δt Between Starts of Repetitions:

- **Frame Spacing** – Distance between frames (i.e. out-of-plane spacing).
- **Scan Speed** – Speed at which transducers move. This setting cannot be controlled by the user.
- **Default ROI Width/Length** – Default width and length of 3D ROIs that can be selected using the **Default** button in the acquisition table in 3D tab. Imaging mode specific.
- **Remember ROI Size and Position** – When checked, the ROI size and position for a given imaging mode will be remembered after switching modes.
- **Step Mode** –
 - **Stepped** – Acquire one or more frames per position. The transducer comes to a complete stop with each step before advancing to the next position.
 - **Frames Per Position** - number of frames to acquire at each step. When used for a static acquisition, multiple frames in the same spatial location will be averaged together during volume fusing (which can improve image signal-to-noise ratio).
 - **Continuous** – Acquire a single frame per position. Transducer sweeps across each pass without stopping.

- **Sequence Duration** – Duration of M-Mode and 2D sequences. Also represents the duration of sequence at each step in a 3D sequence scan.
- **Advanced Timing Settings** –
 - **Δt between stepped positions** – Delay between steps when **Step Mode** is set to “Stepped”.
 - **Δt between start of passes** – Delay between the start of each pass for a multi-pass scan.
 - **Repetitions** – Number of times to repeat a scan. Data from multiple repetitions is averaged together.
 - **Δt between start of repetitions** – Delay between the start of each repetition for a multi-repetition scan.

4.3.2 Post-Scan Settings

Ultrasound Post-Scan Settings

Post-Scan settings control actions that can happen after data acquisition, such as Fusing, which is the process of stitching multiple passes into a single 3D volume.

- **Processing Options** –
 - **Process Now** – Process (3D reconstruction, Phase Gating, SWE processing) scan immediately after acquisition.
 - **Process Later** – Delay processing for a later time. You will be prompted to process any unprocessed volume when exiting SonoEQ. Furthermore, you can process volumes manually by using the tools menu (Tools → Volume Reconstructor, Run Phase Gating, Run SWE Processing).
- **Auto-generate MIP** – write a .PNG file containing a maximum intensity projection (MIP) of the acquired 3D scan
- **Threshold for Breathing Motion** – 3D B-Mode and AA scans can be processed to remove breathing artifacts after acquisition. The threshold for the algorithm (Low, Medium, High) defines how aggressively to remove artifacts with “High” being the most aggressive, which might result in removing data without breath artifacts.
- **Unprocessed Data Storage** –


revvity	SonoEQ Version: 2.0.2		Revision: 22
	Title: SonoEQ User Guide		
	Proprietary & Confidential	Effective Date: 11/22/2023	Page 29 of 108

- **Uncompressed** – Write the unfused and unprocessed data taken with each acquisition without compression.
- **Compressed** – Write the unfused and unprocessed data taken with each acquisition with compression to reduce file size.
- **Remove after processing** – Remove the unfused and unprocessed data after the initial fusing and processing for that data is complete.
- **Fused Volume Z Spacing** – Distance between frames in reconstructed 3D scan.
 - **Same as acquired** – value specified by **Frame Spacing** in Scan-Plan settings.
 - **Custom** – use a different value than that specified by **Frame Spacing** setting. If the reconstruction spacing is set larger than the acquired frame spacing, overlapping frames will be averaged. If the reconstruction spacing is set smaller than the acquired frame spacing, there may be “holes” in the reconstructed image, i.e. black areas where no real data exists to fill the voxel.
- **2D Phase Correction** – enabling this setting processes an acquired sequence to produce a single-cycle cardiac sequence. For more details, see section 2.3.3.



NOTE: 3D sequences are automatically processed with phase correction.

4.3.3 B-Mode Presets



Optimized – Acquire maximum quality images at the cost of longer scan times.

- Suitable for breathing correction
- 0.1 mm spacing (Linear Array), 0.13 mm spacing (Dual Element Wobbler)
- Slower scan speed


Multi-focal Optimized – Acquire images with maximum quality throughout depth.

- 3 focal zones
- Slower scan speed

Fast Scan – Quickly acquire medium quality images.

- No breathing correction
- 0.15 mm spacing (Linear Array), 0.25 mm spacing (Dual Element Wobbler)
- Faster scan speed

4.3.4 Acoustic Angiography Presets




Density – The settings for this preset allow for the capture of all vascular information including resolvable and sub-resolution vasculature.

- Continuous scan
- Multiple repetitions
- Fast scan speed

Morphology – This preset is designed to isolate signal from resolvable vessels ($\geq 100 \mu\text{m}$) for visualization of vascular architecture.


- Stepped scan
- Multiple images per spatial location
- Slow scan speed

4.3.5 Shear Wave Elastography Presets




Optimized – Best parameters for acquiring a high-quality SWE image.

4.3.6 Custom User Presets



Changes to scan settings can be saved into custom presets per imaging mode. Clicking the **Save Preset As** button brings up a dialog where you can name your preset and save it to disk.

To edit a user-defined preset, click the arrow on the side of the custom preset button to access its context menu. From here you can rename, delete, or overwrite the preset with the current scan settings.



While the factory presets do not set all imaging settings, custom presets save all the **Scan Protocol** and **Post-Scan** settings for that mode. Furthermore, a **Shear Wave Elastography** custom preset saves setting for both **Shear Wave Elastography** and **B-Mode** imaging modes.

5 Analyzing Image Data

The following sections describe different ways of analyzing data. In general, analyzing data involves annotating the data (i.e. “segmenting”) to extract parameters of interest (e.g., volume measurement, image intensity, functional metrics, etc.). Datasets with different dimensionality usually require unique tools to perform measurements. These are described in more detail below.

SonoEQ starts in Analysis mode by default, and can be accessed by clicking the **Image Analysis Tab**.



5.1 Loading Image Data

Data can be loaded into SonoEQ using the Load Dialog by clicking the **Load Volume(s)** button.




NOTE: The default save location for all data captured with SonoEQ is:
%UserProfile%\Documents\SonoEQ

The Load Dialog can be populated by adding volumes using the **Select Folder** and **Select File(s)** buttons.

Select Folder – will open a Windows file explorer window where you can point to the location of a data folder. All available scans on that folder and its subfolders will be added to the Load Dialog.

Select File(s) – will open a Windows file explorer window where you can point to a single scan file (MHA, SEQ.MHD, SEW.NRRD). Only that file and any accompanying files, such as the supporting 2D B-Mode acquired during a 2D SWE acquisition, are added to the Load Dialog.


Select Folder and **Select File(s)** can be used multiple times to populate the Load Dialog.



	Volume
1	<input checked="" type="checkbox"/> Study 1 - 2020_04_01 - 1-1 - [13:51:00] - B-Mode - 3D
2	<input checked="" type="checkbox"/> Study 1 - timepoint 1 - 1-1 - [10:17:05] - B-Mode - 2Dt
3	<input checked="" type="checkbox"/> Study 1 - timepoint 1 - 1-1 - [15:05:28] - B-Mode - 3D
4	<input checked="" type="checkbox"/> Study 1 - timepoint 1 - 1-1 - [15:07:13] - B-Mode - 3D
5	<input checked="" type="checkbox"/> Study 1 - timepoint 1 - 1-1 - [08:45:08] - M-Mode - Live
6	<input checked="" type="checkbox"/> Study 1 - timepoint 1 - 1-1 - [08:46:32] - M-Mode - Live
7	<input checked="" type="checkbox"/> Study 1 - timepoint 1 - mouse 1 - [16:17:17] - AA - 3D
8	<input checked="" type="checkbox"/> Study 1 - timepoint 1 - mouse 1 - [11:13:51] - AA - 3D
9	<input checked="" type="checkbox"/> Study 1 - timepoint 1 - mouse 1 - [10:46:38] - B-Mode - 2D

Clear Table File Count: 39 Total File Size: 1418.5 MB OK Cancel

Clicking **OK** will only load the files that are visible and checked in the table. Further tuning can be done by using the Filter text box and checking/unchecking entries in the table. In this example, only the “animal1” 2D entries that are checked will be loaded.



	Volume
1	<input type="checkbox"/> Study 1 - 2020_04_01 - 1-1 - [13:51:00] - B-Mode - 3D
2	<input type="checkbox"/> Study 1 - timepoint 1 - 1-1 - [10:17:05] - B-Mode - 2Dt
3	<input type="checkbox"/> Study 1 - timepoint 1 - 1-1 - [15:05:28] - B-Mode - 3D
4	<input type="checkbox"/> Study 1 - timepoint 1 - 1-1 - [15:07:13] - B-Mode - 3D
5	<input type="checkbox"/> Study 1 - timepoint 1 - 1-1 - [08:45:08] - M-Mode - Live
6	<input type="checkbox"/> Study 1 - timepoint 1 - 1-1 - [08:46:32] - M-Mode - Live
7	<input checked="" type="checkbox"/> Study 1 - timepoint 1 - mouse 1 - [16:17:17] - AA - 3D
8	<input checked="" type="checkbox"/> Study 1 - timepoint 1 - mouse 1 - [11:13:51] - AA - 3D
9	<input checked="" type="checkbox"/> Study 1 - timepoint 1 - mouse 1 - [10:46:38] - B-Mode - 2D

Clear Table File Count: 31 Total File Size: 1118.3 MB OK Cancel


Furthermore, data can be quickly loaded by **Drag-and-Dropping** a file or folder into SonoEQ. Drag-and-dropping a folder will load all available image files in the folder and drag-and-dropping a file will load that file, and any supporting multi-stream images (e.g. B-mode 2Dt from an M-Mode acquisition).

5.1.1 Loading Multi-Modal Data



SonoEQ supports the loading of data acquired with other Revvity In Vivo imaging devices. Loading data from other devices follows the same workflow outlined in Section 5.1.

The following data acquired with the IVIS® Lumina or Spectrum imaging systems and exported with the Living Image® Analysis Software can be viewed in SonoEQ image analysis.

<u>File Type</u>	<u>Extension</u>
Luminescent (2D/3D)	.tiff
CT (3D)	.dcm
Photograph (2D)	.png



Sequences of Luminescent Photographs (2D) can be viewed using the **Playback** widget located in the toolbar. Playing the sequence is controlled via the **Play/Pause** button. Advancing the movie forward or backward can be done by clicking the **Backward/Forward** buttons. The current frame of the movie can be manipulated by adjusting the **Playback Slider**.




The following data acquired with the Quantum GX2 imaging system can be viewed and segmented in SonoEQ image analysis.

File Type

Extension

CT (3D)

.vox






NOTE:

Data acquired with the IVIS® systems must be exported as float-corrected images and be accompanied by a SequenceInfo.txt or ClickInfo.txt file to load. See *Living Image User Manual* for more details on export features.


5.1.2 Viewing Multi-Modal Data Side-By-Side


Data acquired with the Vega can be visualized side-by-side with data acquired with other Revvity In-Vivo imaging devices in SonoEQ. This feature allows for multimodal data comparison of the same mouse.




Vega Ultrasound Data (top) and QuantumGX2 Data (bottom)

To view data from different imaging modalities, first ensure that your two volumes are the selected primary and secondary volumes, then choose the **Side-By-Side** or **Three over three** views from the **Slice View** settings in the toolbar.




To select a non-default slice orientation, hover over the  icon and select the orientation you wish to appear in the slice.




5.2 Viewing/Navigating Image Data

A 3D B-Mode will be used to illustrate how data is displayed and oriented. Upon loading a 3D data set, the Slice Views window in SonoEQ will display a 2x2 grid of windows. This “Four-Up” view shows the three 2D orthogonal views - axial (red), coronal (green), and sagittal (yellow) - along with all three slices simultaneously displayed in the 3D view (blue).




These three orthogonal planes correspond to the Vega imaging system as follows:



NOTE: The ultrasound beam is always coplanar to the Red/Axial plane.

Prior to segmenting the image, it is helpful to quickly scroll through the image slices and orient yourself to common anatomical landmarks. Scrolling can be done by hovering the mouse over a given slice and using the mouse wheel, or by clicking and dragging the slider above each slice as indicated below. Additionally, the image can be zoomed by right-clicking and dragging up or down or panned by clicking and dragging the mouse-wheel button. See Chapter 10 for additional tips.



5.3 Overlaying Multiple Datasets

Vega instruments offer the capability of capturing multi-modal datasets from the same subject (e.g. B-Mode and Acoustic Angiography). These images can be analyzed separately or, if desired, analyzed together by overlaying the captured data. This is accomplished through the **Primary** and **Secondary** selector drop down menus.

Primary

Secondary

- 1) Select two files in the Load File Dialog that represent images captured from the same mouse and same timepoint.

Select Folder

Select File(s)


Filter...

	Volume Paths
1	<input checked="" type="checkbox"/> PKD - Timepoint1 - 60_A - [10:46:59] - AA - 3D
2	<input checked="" type="checkbox"/> PKD - Timepoint1 - 60_A - [10:45:15] - B-Mode - 3D - Corrected: High




NOTE: If the file chosen represents a two-stream dataset (e.g. B-Mode + M-Mode), both datasets will be pulled into the Load File Dialog automatically.

- 2) Select each volume to the Primary and Secondary menus, respectively. Typically, the B-Mode (used for anatomical context) is used as the Primary volume while the advanced imaging mode (e.g. Acoustic Angiography) is used as the Secondary volume.





- 3) By default, the two volumes will be rendered in **Overlay Mode**, with the Secondary volume blended with the Primary volume. Adjust the **Opacity Slider** to set to the opacity of the layers to accentuate either the primary or secondary volume as desired. Below shows an example of Acoustic Angiography image (yellow pixels) overlaid with 50% opacity on top of a B-Mode (gray pixels) positioned over the mouse kidneys.




CAUTION:


When overlaying two volumes, an assumption is made that the animal did not move between scans. If properly sedated, animals should remain in the same position to allow overlaying consecutively acquired data sets. However, if the mouse changes position during imaging (e.g. wakes from anesthesia), the image data will be misaligned in the overlay view.

- 4) Alternatively, Primary and Secondary volumes can be viewed in **Side-By-Side Mode** by clicking the **View Toggle** button next to the Primary volume selector . When viewing volumes in side-by-side, the Slice View for the Secondary volume will be color coded with a


faded color palette and is referred to as the “+” layout (e.g. Red and Red+). Clicking the View Toggle button a second time will return to Overlay View .



NOTE:


Depending on the acquisition parameters, the Primary and Secondary volumes may not span the same spatial extent. For example, a B-Mode may be taken as a widefield scout scan of the whole abdomen, but Acoustic Angiography may be targeted to the kidneys specifically. If this is the case, the Slice views will limit the visible extent to the bounds of the Primary volume. To quickly toggle the positions of the Primary and Secondary volumes (and thereby change the maximum visible extent), the **Switch** button can be utilized .

A pair of 3D static volumes that have overlapping spatial extents support the option to link slice intersections between both volumes. Linking/unlinking slice intersections for a pair of 3D static volumes can be controlled through the link button in the slice controller (visible when mouse hovering over the slice header).



This links the slice offsets of both volume such that adjusting one will adjust the corresponding offset in the other volume (red slice offset in primary volume will change the

offset in the red slice in the secondary volume), which allows for easily viewing the same world location across two datasets such as a B-Mode and Acoustic Angiography volume.



NOTE: Holding shift and moving the mouse cursor over one of the slice viewers will center the other slices on that location. If the volumes are linked, the slices for both volumes are centered on the mouse cursor location.

5.4 Creating Segmentations and Quantifying Data

To quantify various aspects of image data, a “segmentation” must be created. The process of segmenting (also known as contouring) involves delineating structures of interest in 3D space. Due to the wide breadth of applications that are supported by the Revvity instruments, SonoEQ includes a variety of segmentation tools, also known as “segmentation effects”, that can be utilized. The most basic effects mimic a manual painting interface (akin to Adobe Photoshop) but work on 3D matrices of voxels rather than on 2D pixels. Other, more sophisticated effects automate the segmentation process by making assumptions about the geometry of the 3D object that is being painted or utilizing information from the volume itself to simplify the process.

The subsections below are organized such that segmentation effects are introduced first, followed by step-by-step tutorials on how to apply each based on the given datatype and dimensionality (starting with the basic static 3D volume). Note that not all effects are applicable to all datatypes, and SonoEQ will limit the available effects as necessary (e.g. M-mode sequences have different segmentation effects than 3D volumes).


5.4.1 Segmentation Effects Overview

The following segmentation effects are available in SonoEQ. From left to right, they are None, Surface Cut, Draw, Scissors, Paint, Erase, Fill Between Slices, and Smoothing:




Segmentation effects toolbar

- **None** – Deselects the segment effect. Use this option when manipulating the segmentation is not desired (e.g. panning through the data, zooming, 3D render, etc.).
- **Surface Cut** – Uses control points to fill a segment following a convex hull-like algorithm. This effect can be used for quickly segmenting simple spheroidal 3D shapes with little or no concavities, such as a tumor.




- **Point Placement Button** – Toggles the cursor into point placement mode for placing control points. This button must be un-selected to move already placed control points. Right-clicking in a slice view will un-select this button.
- **Delete Point Button** – Deletes the most recently placed control point. Click the dropdown button to show a menu that allows deleting all points in one action.
- **Edit** – Used to re-enter the “Editable” mode of a surface cut segmentation. This button is utilized when a surface cut needs to be modified after clicking the Apply button.
- **Operations** –
 - **Replace (default)** – replaces the old 3D segmentation with the new one.
 - **Fill inside** – adds additional area to the old segmentation.
 - **Erase inside** – deletes overlapping area.
- **Smooth Model** – smooths the resulting segmentation and removes sharp corners. Turning this option OFF will produce a true convex-hull segmentation.
- **Draw** – Creates a segment in a single 2D slice.




- **Left click** – Adds successive points to the drawing region.
- **Left click and drag** – Draws a continuous curve (similar to a lasso tool.)
- **X key** – Delete the last added point.
- **Right click or Enter** – Completes the segmentation and fills the lassoed space.

- **Scissors** – Cut through slices from a single 2D slice.




- **Left click and drag** – Depending on selected operation, draw a continuous curve to cut (similar to a lasso tool) or draw a shape.
- **Esc** – Cancel an unfinished cut.
- **Operation** –
 - **Fill inside (default)** – adds additional area within the lasso.
 - **Fill outside** – adds additional area outside of the lasso.
 - **Erase inside** – deletes overlapping area within the lasso.
 - **Erase outside** – deletes overlapping area outside of the lasso.
- **Slice Cut** –
 - **Unlimited** – Project the cut through the entire volume
 - **Positive** – Only positive side of the slice plane is included in cut region
 - **Negative** – Only negative side of the slice plane is included in cut region
 - **Symmetric** – Cut region is limited to specified thickness around the slice plane

- **Paint** – Paints with a circular brush in either 2D or 3D.




- **Left click** – Paints a single circular region.
 - **Left click and drag** – Paint in a continuous manner.
 - **Shift + mouse wheel** or **+/- keys** – adjust the brush size.
 - **Ctrl + mouse wheel** – Zoom the slice view in/out.
 - **Diameter** – Sets the size of the circular brush region.
 - **Sphere Brush** – Toggles between a 2D and 3D brush. Note that if Sphere Brush is active, segmentation voxels will be filled out of plane and will not be visible in the current slice view until the viewer is scrolled.
 - **Edit in 3D Views** – Allow painting in 3D view. If enabled, click-and-drag in the 3D view paints in the view instead of rotating the view.
 - **Color smudge** – Select segment by sampling the pixel location where the brush stroke starts. If brush stroke starts in an empty area, then the brush erases highlighted region from the selected segment.
- **Erase** – Erases with a circular brush in either 2D or 3D.




- **Left click** – Erases a single circular region.
- **Left click and drag** – Erase in a continuous manner.
- **Shift + mouse wheel** or **+/- keys** – adjust the eraser size.
- **Ctrl + mouse wheel** – Zoom the slice view in/out.
- **Diameter** – Sets the size of the circular eraser region.

- **Sphere Brush** – Toggles between a 2D and 3D eraser. Note that if Sphere Brush is active, segmentation voxels will be erased out of plane.
- **Edit in 3D Views** – Allow erasing in 3D view. If enabled, click-and-drag in the 3D view erases in the view instead of rotating the view.
- **Fill Between Slices** – Interpolates segmentation between slices. This tool is used in conjunction with other effects (e.g. paint, draw), and will fill gaps where empty slices exist. The segmentation will only be expanded if a slice is segmented but none of the direct neighbors are segmented.



- **Preview:**
 - **Initialize** – Initializes the Fill Between algorithm and produces a preview of the filled region in a transparent color in the slice views.
 - **Auto-update** – Toggles whether the Fill Between algorithm should be re-initialized automatically upon modifying the segmentation (e.g. painting a new slice). If enabled, the preview will adjust each time a new slice is added.
- **Display:**
 - **Inputs** – Adjusts the transparency of the input segmentation to be more or less visible based on the position of the slider.
 - **Show 3D** – Toggles whether the segmentation is visible in the 3D view.

- **Smoothing** – Makes segment boundaries smoother by removing extrusions and filling small holes. This effect is applied AFTER an initial segmentation has been drawn using the other effects.




- **Gaussian**
 - **Standard Deviation** – Sets the standard deviation value (σ) of the Gaussian smoothing kernel.
- **Median**
 - **Kernel Size** – Sets the size of the median filter kernel.
- **Opening**
 - **Kernel Size** – Extrusions smaller than this size will be removed.
- **Closing**
 - **Kernel Size** – Sharp corners and holes smaller than this size will be filled.
- **Diameter** – Sets the size of the circular eraser region.
- **Sphere Brush** – Toggles between a 2D and 3D eraser. Note that if Sphere Brush is active, segmentation voxels will be erased out of plane.
- **Edit in 3D Views** – Allow erasing in 3D view. If enabled, click-and-drag in the 3D view erases in the view instead of rotating the view.

The selected segmentation effect is displayed next to the cursor when hovering over slice views to avoid confusion.

5.4.2 Creating a Segmentation with Surface Cut


The following section gives a tutorial on how to draw a segmentation using the “Surface Cut” technique. This approach allows rapid segmenting of 3D shapes with very few clicks. In general, this method works best when features of interest have simple spherical/ellipsoidal geometries (e.g. tumors, kidney, bladder, etc.). As an example, a tumor is segmented below.

- 1) With the feature of interest identified in all three planes, click the **Add** button, and a new row will appear in the “Segmentations” table,



When the **Add** button in the create tab is clicked, the default segmentation effect that is selected is **Surface Cut**. **Surface Cut** allows points to be placed by clicking on the 3D Image slice views. A surface is generated which intersects as many of the placed points as possible.

The goal of the next several steps will be to define the boundaries around the perimeter of the tumor such that a 3D surface (i.e. a segmentation encompassing the entire volume of the tumor) can be produced.



The figure above is a high-level overview of the approach and illustrates the workflow for **Surface Cut**. In this figure, there is no image data shown, but the intention is to extract a feature of interest (e.g. a kidney) and that these points fall at the interface between the tissue of interest, and surrounding tissue.


A) The first three points (1,2, and 3) have been added to the perimeter of the structure of interest. The software then automatically draws a smooth contour encompassing those points.

B-C) As additional points (4 and 5) are added to refine the shape of the surface, the contour will update in real time.


D) If a point is added beyond the existing perimeter bounds, the contour will adapt to include that new point (6, indicated by arrow), disregarding the contribution from the nearby point (4, indicated by red X) because it now lies within the interior of the surface.

The workflow described above is performed using SonoEQ on 3D ultrasound data, but for simplicity was illustrated in 2D in the above figure. The following steps will describe how to perform the Surface Draw procedure on actual 3D ultrasound data.


- The software will automatically enter **Point Placement** mode . You can begin clicking on one of the three image planes to place segment points.




NOTE: Placing points in a single slice view will generate a single-slice segmentation. Placing points in orthogonal slice views (i.e. out-of-plane) will expand the segmentation from 2D to 3D.

- If you made a mistake and need to adjust a point, simply toggle the **Point Placement** mode and click-and-drag points to adjust them.
- To **Delete** a point placement, utilize the trash can button  to delete the most recently placed point. Alternatively, clicking the dropdown arrow on the trashcan button will allow you to **Delete All Points**, and start over.


- 5) Continue placing points around the tissue of interest (e.g. tumor) until the contour lines of the segmentation match the boundary of the object being segmented. The more points that are placed, the more accurate the 3D surface will be. In the following example, 7 points in the Axial view and 2 points in the Coronal view define the tumor region.




- 6) After all points have been placed, click the **Apply** button to fill the region defined by the surface and generate a 3D segmentation object.



- 7) The interpolated segmentation will then appear as a translucent overlay in the slice views and will be rendered in 3D as a translucent surface. Additionally, the volume of the segmentation will be displayed in the Segmentations table.



- 8) After reviewing the segmentation, if you are unsatisfied with its accuracy, simply click the **Edit** button to return into the Surface Cut effect and add/modify/remove points.



Three options determine how an edited segment will appear after clicking Apply:



NOTE:


- **Replace** (default) – replaces the old 3D segmentation with the new one.
- **Fill inside** – adds additional area to the old segmentation.
- **Erase** – deletes overlapping area.

Using the Erase option allows the generation of concave regions in a segmentation.

5.4.3 Creating a Segmentation with Slice Painting


The following section gives a tutorial on how to draw a segmentation using the “Slice Painting” technique. This approach allows more flexibility compared to the Surface Cut approach but is more time intensive. As an example, the liver of a mouse with NAFLD is segmented below.


- 1) With the feature of interest (in this example, assumed to be a liver) identified in all three planes, click the **Add** button.



A new row will appear in the “Segmentations” table,

	Visibility	Opacity	Markup Name	Markup Type	Completed
1		40 %	Markup_1	3D Segmentation	

- 2) Click the **Draw** effect. The software will enter **Point Placement** mode , and you can begin clicking on one of the three image planes to place draw points. Clicking and dragging will create a continuous line that follows the mouse cursor, or alternatively individual clicks will connect with a straight line.
- 3) To **Delete** a point placement, press the “x” key on the keyboard.
- 4) Continue placing points around the tissue of interest (e.g. liver) until the contour lines of the segmentation match the boundary of the object being segmented. Once complete, either click the **Right mouse button** or press “Enter” to fill the region.



Original image

Outlined 2D region

Filled 2D region

- 5) Advance the slice view using the mouse wheel or the slice slider and repeat the process, drawing sequent 2D outlines to cover the length of the feature of interest. If desired, each individual 2D frame can be segmented manually following this process, however this is typically very time consuming. Instead, it is advisable to only segment a subset of the slices and use the **Fill Between Slices** and **Smoothing** effects as described below.
- 6) After drawing several 2D slices and covering the extent of the feature of interest, navigate to the **Fill Between Slices** effect and click the **Initialize** button. Once initialized, the interpolation algorithm will be run automatically filling the empty space between the regions. New 2D segmentations can be added as desired. Once finished, click the **Apply** button to complete the segmentation.



1 slice painted

5 slices painted


Filled slices




NOTE:

Filling between slices will only interpolate and does not do any extrapolation. This means that one slice at the top and bottom of the feature of interest must be drawn, otherwise it will appear as a truncated segmentation.


- 7) (Optional) Sometimes after the slices are filled, there are sharp transitions that may be undesirable for visualization purposes. The 3D filled segmentation can be smoothed using the **Smoothing** effect.



Original segmentation after fill




Smoothed segmentation



3D rendering

5.4.4 Quantifying a Segmentation

To quantify the segmentations and export results to disk, click on the **Quantify** tab within the Control Panel.



Measurements are broken into “Surface statistics” and “Image statistics”. Surface statistics represent gross parameters of the segmentation object such as volume and surface area.

Surface statistics				
	Voxels	Volume (mm ³)	Volume (cc)	Surface Area (mm ²)
Markup_1	5089	1.963	0.002	38.136

Image statistics apply the segmentation as a mask to the underlying volume and compute parameters such as mean pixel intensity. Image statistics include an optional **Threshold Slider** that will remove pixels that fall outside of the selected upper and lower bounds. By using the threshold, a Pixel Positivity metric can be computed, which represents the number of pixels in the segmentation that fall above the threshold divided by the total number of pixels.

▼ Image statistics (Primary)


	Positivity (%)	Mean (count)	Median (count)	Standard Deviation (count)
Markup_1	100	81.880	82	28.599

☒ Apply threshold to statistics

Threshold:
Manual

0
255

In the case of an “overlapping” dataset (where Primary and Secondary volumes are selected), the image statistics will be computed for each using the same segmentation object.


You can save your work at any point by clicking on the **Save** button  located in the toolbar. To resume your work in another SonoEQ Analysis session, simply load in the volume(s), and the last saved state will be loaded depending on the elements present in the last saved session. Frequent saving is recommended when working with a large batch of volumes.


5.4.5 Analyzing Data with Digital Calipers


Digital calipers offer a quick solution to measuring linear distances in the 3D datasets captured by SonoEQ. Calipers can be placed in either Acquisition or Analysis mode.

- 1) To start placing a caliper, click the **Add Caliper** button located in the top right corner of the SonoEQ interface



- 2) The software will enter **Caliper Placement Mode**  and you can begin placing the start and end points of the caliper by clicking within a slice view or the 3D view where you would like to measure. The first click places one end of the caliper, and the second places the other. After the first click, a yellow ruler will appear with the distance of the caliper labeled based on where their mouse is currently hovering. After a second click, the second end point of the caliper is placed. Click the Add Caliper button again to place multiple calipers.



- 3) If you made a mistake and need to adjust the caliper, move the mouse over one of the crosshair endpoints and click-and-drag to adjust. The cursor will turn into a hand icon  to indicate the point is being edited when dragging.



CAUTION:

Volume measurements derived from calipers tend to be less accurate compared to volume estimates derived from 3D segmentations. This is because the equations that convert calipers measurements into a volume typically assume an underlying geometry (e.g. an ellipsoid) which may not accurately reflect the tissue of interest.


- 4) To delete an individual caliper, click the **Delete Calipers** button in the upper right corner of the SonoEQ interface. A message box will appear.



- 5) Caliper visibility can be toggled on and off with the **Show/Hide Caliper** button located in the upper right corner of the SonoEQ interface.




- 6) To quickly locate and view calipers placed in different positions within the scene, use the **Select a Caliper** drop-down menu, and the slices will update to the position of the selected caliper.



- 7) To save calipers to disk, click on the Save button located in the top left part of the interface. Upon loading the volume in a future SonoEQ session, any calipers that were saved will be automatically loaded with their corresponding volume.



- 8) Lastly, to export caliper measurements to disk, click on **File → Export Statistics from Scene → Caliper**. A save file dialog will appear and allow you to save the caliper measurements as a CSV file.



5.4.6 Analyzing B-Mode Sequences


Segmenting B-mode sequences proceeds in much the same way as segmenting a static image. The primary difference is that segmentations are linked to a given frame in the movie sequence.

The following describes the procedure for analyzing the volume of the left ventricle (LV) to assess functional metrics such as Ejection Fraction (EF) during cardiac imaging.

- 1) Once in Analysis mode, the first step is to load one or more images for review. Clicking the **Load Volume(s)** button will bring up a file dialog that allows you to select either individual files or entire folders of data.



- 2) Upon loading a sequence data set, the Slice Views window in SonoEQ will display the traditional four-up view showing axial, sagittal, coronal, and 3D views.
- 3) (Optional) If the ultrasound planes do not match up with the long and short axis of the heart, the ortho-slices can be rotated and re-sliced for a better match. Press **Ctrl + Alt + Left-click-and-drag** in any slice to rotate the intersections. This step is useful to better identify anatomical landmarks, such as the papillary muscles, and makes segmentation easier to perform. An image of a re-sliced volume is shown below.



- 4) Prior to segmenting the image, it is helpful to play the sequence movie and identify frames depicting peak systole and end diastole. This can be done with the **Playback** widget that is in the toolbar. Playing the movie is controlled via the **Play/Pause** button. Advancing the movie forward or backward can be done by clicking the **Backward/Forward** buttons. The current frame of the movie can be manipulated by adjusting the **Playback Slider**.



- 5) Click the **Add** button to create a new “3D Segmentation Sequence” markup.



A new row will appear in the Markup table, that includes two sub-rows corresponding to the four traces that need to be drawn to fully quantify cardiac parameters.

The goal of the next several steps will be to define the boundaries of the epi- and endocardial surfaces that will be used to compute LV function.

	Visibility	Opacity	Markup Name	Markup Type	Completed
▼ <input type="radio"/>	<input checked="" type="radio"/>		Markup_1	3D Segmentation Sequence	
1 <input type="radio"/>		40 %	3 <input checked="" type="checkbox"/> Endocardial Sys/Dia	Frame 33	2
<input type="radio"/>		40 %	... Endocardial Sys/Dia		
<input checked="" type="radio"/>		40 %	... Epicardial Sys/Dia		
			... Epicardial Sys/Dia		

The radio buttons to the left of each sub-row (1) allow you to select the appropriate cardiac surface (e.g. endocardial or epicardial). The drop-down menu on the right (2) indicates which frame in the movie sequence the segmentation is assigned to. Upon placing two segmentations in different movie frames, the software will automatically distinguish systolic and diastolic phases (3) by comparing the sizes of the regions.

The following four segmentations need to be completed for a full assessment of the LV:

- Endocardial volume at peak systole (*LVVs*)
- Endocardial volume at end diastole (*LVVd*)
- Epicardial volume at peak systole (*EpiLVVs*)
- Epicardial volume at end diastole (*EpiLVVd*)

Ejection fraction (EF), stroke volume (SV), cardiac output (CO), and LV mass (LVM) can then be calculated using the following equations given a known heart rate (HR) and assumed cardiac tissue density of 1.05 mg/μL:


$$EF (\%) = \left(\frac{LVVd - LVVs}{LVVd} \right) \times 100$$

$$SV = LVVd - LVVs$$


$$CO = SV \times HR$$

$$LVM = 1.05 \times (EpiLVVd - LVVd)$$


- 6) Using the Playback Slider, select the desired frame in the movie file (either peak systole or end diastole), and start by segmenting the **endocardial** border of the LV up to level of the mitral valve using any of the Segmentation Effects described previously.




If desired, toggle the radio button, and segment the **epicardial** border of the LV at the same frame position. Epicardial segmentations are only necessary if LV mass is being quantified. The epicardial border is not required for functional outputs such as Ejection Fraction (EF) or Cardiac Output (CO).










- 7) Next, use the Playback Slider to select the opposite cardiac phase (in this case peak systole, as end diastole was segmented first above). Repeat the process of segmenting both the endocardial and epicardial borders, using the appropriate radio button selections.



NOTE:

Point placement mode can be turned off quickly by right-clicking, pressing spacebar or by clicking the **Point Placement** mode button .

- 8) Once all segmentations are placed, green checkmarks will appear next to the **Endocardial** and **Epicardial** labels in the table, and a quick preview measurement (e.g. EF%) will be displayed in the Segment table. Quick toggling between systole and diastole can be achieved by selecting the corresponding measurement frame from the drop-down menu.

 Add		 Add (auto)		 Delete	
	Visibility	Opacity	Markup Name	Markup Type	Completed
▼ 			Markup_1	3D Segmentation Sequence	✓ EF (%): 49.947
		40 %		✓ Endocardial Sys ✓ Endocardial Dia	Frame 80 ▼
		40 %		✓ Epicardial Sys ✓ Epicardial Dia	Frame 33 Frame 80

- 9) To quantify all segmentations and export results to disk, click on the **Quantify** tab. LV volume in systole and diastole, ejection fraction (EF), stroke volume (SV), etc. will be computed automatically and displayed within the table.

Create

Quantify

Utilities

Display Options

Info

▼ Cardiac measurements

	LV Vol;d (uL)	LV Vol;s (uL)	LV EpiVol;d (uL)	LV EpiVol;s (uL)	Frame;d	Frame;s
Markup_1	61.522	32.152	82.732	70.832	33	81


◀

▶

▼ Cardiac calculations

	EF (%)	SV (uL)	HR (BPM)	CO (mL/min)	LV Mass (mg)
Markup_1	47.740	29.370	361.000	10.603	86.868

Export All Statistics

- 10) You can save your work at any point by clicking on the save button  located above the Slice Views. To resume your work in another SonoEQ Analysis session, simply load in the volume(s), and the last saved state will be loaded depending on the elements present in the last saved session.


5.4.7 Analyzing M-Mode Sequences

Analyzing M-Mode data of the left ventricle (LV) involves segmenting/tracing the cardiac walls through the cardiac cycle. As with 3D Segmentations, this process is completed in the **Image Analysis** tab.


- 1) Once in Analysis mode, the first step is to load one or more images for review. Clicking the **Load Volume(s)** button will bring up a file dialog that allows you to select either individual files or entire folders of data.



- 2) Upon loading a data set, the Slice Views window in SonoEQ will display a dual view showing the M-Mode trace with a time scale bar below and corresponding B-Mode sequence. The green triangle overlaid on the left side of the image represents the transducer focus at the time of acquisition and is displayed for reference.



- 3) Prior to segmenting the image, it is helpful to quickly scroll through the M-Mode movie frames and select a measurement region not corrupted by breathing. This can be done with the **Playback** widget that is in the toolbar.





NOTE: B-Mode and M-Mode are both saved as “Sequence” files. When two sequences are loaded into SonoEQ, two Playback widgets will appear in the toolbar that control the B-Mode and M-Mode respectively.

- 4) After selecting the appropriate frame in the M-Mode movie file, click the **Add** button to create a new “Cardiac Curve” markup.





A new row will appear in the Markup table, that includes four sub-rows corresponding to the four traces that need to be drawn to fully quantify cardiac parameters.


	Visibility	Opacity	Markup Name	Markup Type	Completed
▼ ●		100 %	Markup_1	CARDIAC CURVES	
●				Posterior Epicardial	
○				Posterior Endocardial	
○				Anterior Endocardial	
○				Anterior Epicardial	

The goal of the next several steps will be to define the boundaries of the epi- and endocardial surfaces that will be used to compute LV function. The radio buttons to the left of each sub-row allow you to select the appropriate wall as follows, where anterior and posterior are defined relative to the ultrasound transducer.

LV Short-Axis View



← Posterior Epicardial
 ← Posterior Endocardial
 ← Anterior Endocardial
 ← Anterior Epicardial

- 5) The software will bring up the **Open Curve Draw** module and automatically enter **Point Placement** mode . Begin by clicking on the M-Mode image to place control points. For

best results place at least one control point at each point in peak systole and peak diastole for at least 4-5 cycles (the more cycles drawn, the more accurate the final measurement).

Open curve draw

Place points in successive order and press apply to generate the Cardiac Curve.

Point Placement:


Edit


Cancel

Apply

- 6) If you made a mistake and need to adjust a point, simply toggle the **Point Placement** mode off and click-and-drag points to adjust them.




NOTE: Point placement mode can be turned off quickly by right-clicking, pressing spacebar or by clicking the **Point Placement** mode button .


- 7) To **Delete** a point placement, utilize the trash can button  to delete the most recently placed point. Note that deleting all points will result in clearing the edits made to the curve.
- 8) After all points have been placed, click the **Apply** button or press the “a” key to accept the line. The software will automatically step to the next cardiac curve and enter the Point Placement mode.

Apply

Continue until all four curves have been drawn to compute all cardiac parameters. After completing the two endocardial curves, the software will automatically detect peak systole and peak diastole and indicate the selections by showing red and green dotted vertical lines.



Editing Mode



Complete



NOTE: If only interested in calculating cardiac output parameters (e.g. ejection fraction, fractional shortening), just the endocardial borders need to be segmented. The epicardial borders are only used in calculating LV mass.

- 9) To quantify the segmentations and export results to disk, click on the **Quantify** tab within the Control Panel. LV inner diameter (LVID), wall thicknesses (LVPW, LVAW), ejection fraction (EF), etc. will be computed automatically and displayed within the table.

Create	Quantify	Utilities	Display Options	Info
--------	----------	-----------	-----------------	------

Cardiac measurements						
	LVID;d (mm)	LVID;s (mm)	LVPW;d (mm)	LVPW;s (mm)	LVAW;d (mm)	LVAW;s (mm)
Markup_1	2.755	1.354	1.323	1.969	1.310	1.800

Cardiac calculations								
	EF (%)	FS (%)	SV (uL)	HR (BPM)	CO (mL/min)	LV Vol;d (uL)	LV Vol;s (uL)	LV Mass (mg)
Markup_1	83.699	50.853	23.766	498.694	11.852	28.394	4.629	140.997

The “Cardiac Calculations” are computed from the M-Mode measurements using the classic Teichholz equations that assume an ellipsoidal geometry for the left ventricle (LV) [1]:

$$LVS = \frac{7 LVIDs^3}{(2.4 + LVIDs)}, \quad LVD = \frac{7 LVIDd^3}{(2.4 + LVIDd)}$$

$$EF (\%) = \left(\frac{LVD - LVS}{LVD} \right) \times 100, \quad FS (\%) = \left(\frac{LVIDd - LVIDs}{LVIDd} \right) \times 100$$


$$SV = LVD - LVS, \quad CO = SV \times HR$$

$$LVM = 1.05 \times [(LVIDd + LVPWd + LVAWd)^3 - (LVIDd)^3]$$

$$Teicholz \text{ Correction Factor} = 7 / (2.4 + D)$$

where LVS and LVD = LV end systolic/diastolic volume, LVID;s and LVID;d = LV end systolic/diastolic inner diameter, LVPW;s and LVPW;d = LV end systolic/diastolic proximal wall diameter, LVAW;s and LVAW;d = LV end systolic/diastolic anterior wall diameter, SV = stroke volume, CO = cardiac output, HR = heart rate, LVM = LV mass. A correction factor is also applied to the mass LV calculation (LVM_{corr}) where D equals length of the minor axis.

[1] J. Stypmann, et al., “Echocardiographic assessment of global left ventricular function in mice,” *Lab. Anim.*, vol. 43, no. 2, pp. 127–37, Apr. 2009.

- 10) You can save your work at any point by clicking on the save button  located above the Slice Views. To resume your work in another SonoEQ Analysis session, simply load in the volume(s), and the last saved state will be loaded depending on the elements present in the last saved session.


5.4.8 AI-Assisted Automated M-Mode Analysis

While cardiac wall segmentation can be done by hand, SonoEQ can automatically estimate placement of each trace in M-Mode images using artificial intelligence.


- 1) Once in Analysis mode, click the **Load Volume(s)** button to bring up a file dialog that allows you to select either individual files or entire folders of data.





- 2) Upon loading a data set, the Slice Views window in SonoEQ will display a dual view showing the M-Mode trace with a time scale bar below and corresponding B-Mode sequence. The green triangle overlaid on the left side of the image represents the transducer focus at the time of acquisition and is displayed for reference.



- 3) Prior to segmenting the image, use the **Playback** widget in the toolbar to scroll to the desired frame in the sequence.



- 4) After selecting the desired frame, click the **Add (auto)** button to generate the estimated cardiac wall traces. The segmentations will cover roughly the whole frame and the process will take a few seconds.



- 5) The estimated traces will appear in the Markup table and on the image like manually segmented traces. As such, the same functionality that applies to manually segmented traces applies to the estimated traces; the **Open Curve Draw** module can be used to edit the estimated traces as described in the previous section “**Analyzing M-Mode Sequences.**”

Add		Add (auto)		Delete	
	Visibility	Opacity	Markup Name	Markup Type	Completed
▼ <input checked="" type="radio"/>		100 %	Markup_1	CARDIAC CURVES	✓ EF (%): 80.872
<input type="radio"/>				Posterior Epicardial	
<input type="radio"/>				Posterior Endocardial	
<input type="radio"/>				Anterior Endocardial	
<input checked="" type="radio"/>				Anterior Epicardial	


- 6) To include only a specific range in the frame for calculations, the **Crop Model** widget can be used to crop a section of the generated curves as a new markup. This can be useful for ignoring sections of breathing artifacts. The widget appears below the Markup table once both endocardial curves have been placed.

Generate a new markup by cropping an existing markup.


Create ROI

Apply ROI Crop

- 7) Click the **Create ROI** button to place a region of interest on the image. The ROI can be dragged and adjusted by clicking and dragging the points on the generated box. Once the ROI specifies the desired region, click the **Apply ROI Crop** button. This will generate a new markup group in the Markup table for the cropped curves.



Full trace with Crop object overlaid




Cropped M-Mode


<input checked="" type="radio"/>	<input checked="" type="checkbox"/>	100 %	Markup_2	CARDIAC CURVES	✓ EF (%): 79.502
<input type="radio"/>				Posterior Epicardial	
<input type="radio"/>				Posterior Endocardial	
<input type="radio"/>				Anterior Endocardial	
<input checked="" type="radio"/>				Anterior Epicardial	

5.4.9 Analyzing SWE Images


Quantifying tissue stiffness can be accomplished by segmenting the SWE volume in the **Create** tab. The volume can be segmented manually using techniques described previously by pressing the **Add** button and using the Segmentation Effect Tools, or automatically by pressing the **Add (auto)** button, which selects the entire SWE volume.




In addition to the SWE volume, SWE acquisitions compute a confidence map with the same dimensions, in which each pixel value represents the confidence in the shear wave velocity measurement at that location in the volume. The confidence floor can be adjusted in the **Quantify** tab to threshold any pixels with a confidence lower than the floor.



SWE stats are computed only on pixels not excluded by the confidence floor or volume threshold and can be found in the SWE statistics (Primary) table (above). The effect of adjusting the confidence floor on the volume can be seen here:




SWE images and statistics can be represented in three different units; **Shear Wave Velocity (m/s)**, **Shear Modulus (kPa)**, and **Youngs Modulus (kPa)**. You can select your preferred representation in Settings → Analysis → Quantification Settings → SWE Representation.



Each SWE representation has a default color scalar bar range. You can customize the SWE color scalar bar range in Settings → Analysis → Utilities and Display Options Settings → Default SWE Color Scalar Bar.

<u>SWE Representation</u>	<u>Default Color Scalar Bar Range</u>
Shear Wave Velocity (m/s)	0.0 - 2.9
Shear Modulus (kPa)	0.0 - 8.3
Youngs Modulus (kPa)	0.0 - 25.0



Default Range

Custom Range

Recall that the following equations are utilized by SonoEQ to convert shear wave velocity (c_s) into either shear (G) or Young's (E) modulus, assuming a tissue density (ρ) equal to 1 g/cm³ and a Poisson's ratio (ν) equal to 0.5:



NOTE:

$$G = \rho c_s^2$$

$$E = 2G(1 + \nu)$$

$$E = 3G = 3\rho c_s^2$$

5.4.10 Analysis of Secondary Volumes

Whenever a 3D segmentation, a 3D sequence segmentation, or a SWE segmentation is added to a primary volume, Image Statistics are calculated for the secondary volume, if one is selected. For example, if a SWE volume is loaded, the corresponding 2D B-Mode will be automatically selected as the secondary volume and adding a markup to the primary volume will compute the SWE and image statistics for the primary and secondary volumes, respectively.

Primary StudyName - Timepoint1 - Animal1 - [10:38:16] - SWE - 2D

Secondary StudyName - Timepoint1 - Animal1 - [10:38:16] - B-Mode - 2D

Create Quantify Utilities Display Options Info

SWE statistics (Primary)

Confidence floor: 0.01

	Positivity (%)	Mean (m/s)	Median (m/s)	Standard Deviation (m/s)	M
Manual	100	0.133	0.001	0.497	
Auto	100	0.300	0.001	0.677	

☒ Apply threshold to statistics

Threshold: Manual

0.00 4.00

Image statistics (Secondary)

	Positivity (%)	Mean (count)	Median (count)	Standard Deviation (count)
Manual	100	114.700	113	22.556
Auto	100	90.650	91	20.361

☒ Apply threshold to statistics

Threshold: Manual

0 250


5.5 Additional Utilities

5.5.1 M-Mode 360


The following describes the workflow for the **M-Mode 360** module, which allows drawing an arbitrarily placed M-Mode line (also known as “Anatomical M-Mode” in the literature) within a 3D cardiac sequence dataset.





- 1) After loading a 3D sequence, click the **Utilities** tab to access the M-Mode 360 tool.



- 2) (Optional) If the ultrasound planes do not match up with the long and short axis of the heart, the ortho-slices can be rotated for a better match. Press **Ctrl + Alt + Left-click-and-drag** in any slice to rotate the intersections. This step is usually recommended so that the M-Mode 360 line can be placed in a single 2D view, rather than having to place it in the 3D view.



- 3) When satisfied with the 3D orientation of the heart, enter **Point Placement** mode  by clicking the respective button in the left panel. You can then click on any of the slice intersections to place the first control point and click a second time to place the second point that represents the desired M-Mode 360 line. After the first click, the M-Mode view will update dynamically with the mouse movement and the green line will be updated to indicate where the data is coming from.



- 4) After completing the M-Mode line placement, either click the **Save M-Mode** or **Save and Analyze M-Mode** buttons. The **Save M-Mode** button will save the M-Mode to disk while the

Save and Analyze M-Mode button will save the M-Mode to disk and automatically launch the M-Mode analysis interface as described above.



NOTE:

The M-Mode image that is produced by **M-Mode 360** is synthetically produced by repeating the single cycle data from the Phase Gated dataset to fill a standard 2 sec buffer. Phase gated 3D cardiac sequences represent one single heart cycle at a very high synthetic frame rate. The repetitions make segmenting the M-Mode easier with the **Open Curve Tool**.

5.5.2 Volume Projection

In scientific visualization, volume projection is a method to display 3D data in a 2D image. A 2D pixel is generated by performing a mathematical operation (e.g. mean, median, maximum) along parallel rays through the 3D volume traced from the viewpoint to the plane of projection. In SonoEQ, volume projection is particularly useful to create images that contain more information than just a simple slice through the data.




Image of a maximum intensity projection in the coronal plane of acoustic angiography data

▼ Volume Projection

Projection Type:

Maximum

Axis of Projection:

Coronal

Fit ROI to Volume Bounds

Display ROI


View Projection

Save Projection

- **Projection type** – Controls the mathematical operation that dictates the projection. Default is maximum intensity.
- **Axis of Projection** – Decide which Axis the projection should be done to. Default axis is Coronal.
- **Fit ROI to Volume Bounds** – Resets the ROI to encompass the full volume.
- **Display ROI** – Toggles the visibility of the ROI.
- **View Projection** – Show/hide the projection. When the projection is shown, an ROI object appears in each slice. By using the ROI interaction handles, you can control the region being projected. The slice containing the axis of projection will be automatically updated when updating the ROI bounds.
- **Save Projection** – Save a copy of the current projection image to disk.

5.5.3 Data Probe

The data probe utility displays information about the volume data at the location of the mouse pointer in the Slice views. Simply hover the mouse over any pixel in the slice views, and the data probe will update in real time.




The Data probe shows several pieces of information:

- **The first line displays** –
 - Color and name of the Slice view that the mouse pointer is currently over.
 - The patient-space location corresponding to the mouse pointer in Right/Anterior/Superior (RAS) space.
 - The orientation of the slice view.
 - The spacing between slices for the slice view.
- **The lower block of three lines displays** –
 - L, F, or B for Label, Foreground, and Background, respectively.
 - The name of the volume (or None if nothing is selected).
 - The IJK coordinates corresponding to the mouse pointer.
 - The value of the volume at that coordinate (a single value for scalar volumes or the number or components for non-scalar volumes).

5.6 Display Options

SonoEQ offers several tools to help with visualization of 3D image data, including filters, colormaps, and volume rendering. These can be found under the **Display Options** tab when viewing files for analysis.




5.6.1 Color Scalar Bar Display






- **Enable Color Scalar Bar** – enabling this checkbox adds scale bar for the current color map.
- **Color Scale Bar Volume Selection** – drop-down allows in which volume to add the color scale bar.

5.6.2 Volume Display

The Volume Display (Primary) and Volume Display (Secondary) drop-down menus allow you to change how SonoEQ displays image data in the 2D slice views for the primary and secondary volumes, respectively.



- **Lookup Table** – Represents the underlying colormap used to display the data. B-Mode data is defaulted to use the Grey colormap, while acoustic angiography data is defaulted to use the Yellow colormap. Defaults can be changed in the **Settings** dialog.
- **Interpolate** – When checked, slice views will display linearly interpolated slices through input volumes. Unchecked indicates nearest neighbor resampling.
- **Window/Level** – When  the window/level settings for the active volume are locked and can't be changed with either the min/max slider or the cursor tool . When  window/level settings for the active volume can be changed.
- **Window/Level Controls** – Double slider with text input to define the range of input volume data that should be mapped to the display grayscale. Auto window level tries to estimate the intensity range of the foreground image data. On mouse over, a popup slides down to add support for large dynamic range by giving control over the range of the window level double slider.
- **Threshold** – Controls the range of the image that should be considered transparent when used in the foreground layer of the slice display. Same parameters also control transparency of slice models displayed in the 3D viewers.
- **Histogram** – Shows the number of pixels (y axis) vs the image intensity (x axis) over a background of the current window/level and threshold mapping.





NOTE:



Ultrasound data is saved as unsigned 8-bit integers with values ranging from 0-255. These are unitless values referred to as “counts”. Scaling the min and max of the window/level control between these two values represents the full dynamic range of the images.


5.6.3 Volume Rendering

The Volume Rendering drop-down creates a 3D rendering of the image data. 3D renderings are shown in the **3D Viewport** which has a blue accent.




3D rendering of a B-Mode volume with a feature segmented in blue

- **Crop** – Simple controls for the cropping box (ROI). More controls are available in the "Advanced..." section. Enable/Disable cropping of the volume with the checkbox. Show/Hide the cropping box using the Display ROI button . Reset the box ROI to the volume's bounds using the Fit to Volume button .
- **Rendering** – Select a volume rendering method. A default method can be set in the application settings Volume Rendering panel.
 - **VTK CPU Ray Casting**: Available on all platforms, the volume rendering is entirely realized on the CPU.
 - **VTK GPU Ray Casting**: Available on all platforms with a dedicated graphics card, fastest method to volume render (default).
- **Techniques** – Advanced properties of the current volume rendering method.




- **GPU Memory Size**: Amount of memory to allocate on the GPU for volume rendering. By default, all the memory of the GPU is allocated. You can change the default behavior in the application settings Volume Rendering panel.
- **Quality** – Not supported yet.
- **Interactive speed** – Ensure the given frame per second (FPS) is enforced in the views during interaction. The higher the FPS, the lower the resolution of the volume rendering
- **Volume Properties** – Advanced views of the transfer functions.



- **Synchronize with Volumes module –**
 - Click** – Apply once the properties (window/level, threshold) of the Volumes module to the Volume Rendering module.
 - Toggled** – By clicking on the checkbox, you can toggle the button. When toggled, any modification occurring in the Volumes module is instantaneously applied to the volume rendering.
- **Scalar Opacity Mapping** – Opacity transfer function. Threshold mode: controls the transfer function using range sliders in addition to control points.
- **Scalar Color Mapping** – Color transfer function.
- **Gradient Opacity** – Gradient opacity transfer function. This controls the opacity according to how large a density gradient next to the voxel is.
- **Mapping controls –**
 - **Left button click** – Set current point or create a new point if no point is under the mouse.
 - **Left button move** – Move the current or selected points if any.
 - **Right button click** – Select/unselect point. Selected points can be moved at once.
 - **Right button move** – Define an area to select points.
 - **Middle button click** – Delete point under the mouse cursor.
 - **Right/Left arrow keys** – Change of current point.
 - **Delete key** – Delete the point and set the next one as current.
 - **Backspace key** – Delete the point and set the previous one as current.


- **ESC key** – Unselect all points.
 - **Interpolation** – Linear or nearest neighbor interpolation.
 - **Shade** – Enable/Disable shading. Shading uses light and material properties.
 - **Material** – Material properties of the volume to compute shading effect.
- **ROI** – More controls for the cropping box.



- **Display ROI** – Show/Hide the bounds of the ROI box.
- **Interactive mode** – Control whether the cropping box is instantaneously updated when dragging the sliders or only when the mouse button is released

5.6.4 Mean Filter

The mean filter is a tool that can be used to reduce the speckle/graininess of the ultrasound image and improve contrast between organs. This mode works best when no additional speckle reduction techniques are applied (e.g. spatial compounding).



▼ Mean Filter

Radius:


Apply

View Filtered Volume

- **Radius** – Controls the kernel size of the filter in units of mm. The larger the radius, the more smoothing and the blurrier the image will become. Kernel sizes between 0.2 mm and 0.4 mm are ideal.
- **Apply** – Executes the filtering operation. The time to compute the filtered image is proportional to the size of the underlying 3D volume and the chosen radius. The larger the data and radius, the longer the computation time.
- **View Filtered Volume** – Show/hide the filtered volume.

5.7 Viewing Study Info

The Info tab displays metadata saved with the image file that is being viewed. Metadata includes information you entered during image acquisition (e.g. study and animal information) and system parameters and settings (e.g. transducer type, focal depths, etc.).



Filter:

▼ Volume Information (Primary)

Attribute	Value
Animal ID	Animal1
Imaging Mode	B-Mode
Study Name	StudyName
Timepoint	Timepoint1
Timestamp	2019_10_22-11_06_49
Selected Transducer	LA1
Focal Depth 0 (mm)	7.0

▼ Volume Information (Secondary)

Attribute	Value
Animal ID	Animal1
Imaging Mode	B-Mode
Study Name	StudyName
Timepoint	Timepoint1
Timestamp	2020_05_08-10_38_16
Selected Transducer	LA1
Focal Depth 0 (mm)	7.0


The displayed information can be filtered by the following subheadings. Multiple filters can be applied by clicking several subheadings (e.g. Study & Animal).

- **Study** –entered info about the study including study name and timepoint.
- **Animal** –entered info about the animal including strain, model, etc.
- **Transducer** – Settings selected during image capture including transducer type, frequency, focal depths, compounding, etc.
- **Scan plan** – Information about the scanned region.
- **Other** – Software versions and other miscellaneous information.

Certain metadata can also be edited from the table. Metadata Value fields with a white background can be double clicked to begin editing. Edited values will show a (*) by the attribute name indicating that it is unsaved. Edited values are not saved to the image file until the Save button is clicked and the scene is saved.

5.8 Analysis Settings

Analysis settings and defaults can be accessed by clicking **Edit → Settings** from the top menu or opening the **Settings Tray**. Analysis settings are broken down into three subcategories including Markup Settings, Quantification Settings, and Image Processing settings.




- **Create Markup Settings** – Settings related to the segmentation wizard (**Create** tab).
 - **Root name** – Sets the default name for new segmentations. SonoEQ will automatically append a numerical entry (1, 2, etc.) to the given root name when creating new segmentations.
 - **Default colors** – Sets the color pattern order for segmentations. Clicking one of the colored buttons allows for customization of the color.
 - **3D Segmentation Default Opacity** – Sets the default transparency level for the 3D segmentation objects
 - **Cardiac Curves Default Opacity** – Sets the default transparency level for segmentation related to cardiac analysis (e.g. M-Mode curves)
 - **Point Glyph Scale** – Changes the size of the control points when placing new segmentations
 - **Line Width** – Changes the width of the line connecting the control points when placing a new segmentation
 - **Link Opacity of Preview and Applied Markup** – When checked, the opacity settings apply to both the 3D segmentation object, as well as the defined control points. When unchecked, only applies to 3D models and points are set to 100% opacity.

- **Preferred 3D Segmentation Sequence Effect** – Determines which segmentation effect is selected by default when creating a 3D sequence segmentation.
- **Preferred 3D Segmentation Effect** – Determines which segmentation effect is selected by default when creating a 3D segmentation.
- **Preferred 2D Segmentation Effect** – Determines which segmentation effect is selected by default when creating a 2D segmentation.
- **Preferred SWE Segmentation Effect** – Determines which segmentation effect is selected by default when creating a SWE segmentation.

The screenshot shows the 'Analysis' menu on the left and the 'Create Markup Settings' dialog box on the right. The 'Analysis' menu is open, showing options like B-Mode, Acoustic Angiography, Shear Wave Elastography, and Analysis. The 'Create Markup Settings' dialog box is open, showing the 'Quantification Settings' tab. The 'Quantification Settings' tab has several dropdown menus for 'Shown in 3D Segmentation Table', 'Shown in 2D Segmentation Table', 'Shown in Cardiac Curves Table', 'Shown in 3D Segmentation Sequence Table', 'Shown in SWE Segmentation Table', and 'SWE Representation'. The 'Default SWE Confidence Floor' is set to 0.75.

- **Quantification Settings** – Settings related to the quantification table (**Quantify** tab).
 - **Shown in 3D Segmentation Table** – Determines which quantification value is previewed in the **Create** table for 3D segmentations.
 - **Shown in 2D Segmentation Table** – Determines which quantification value is previewed in the **Create** table for 2D segmentations.
 - **Shown in Cardiac Curves Table** – Determines which quantification value is previewed in the **Create** table for cardiac M-Mode segmentations.
 - **Shown in 3D Segmentation Sequence Table** – Determines which quantification value is previewed in the **Create** table for cardiac 4D segmentations.
 - **Shown in SWE Segmentation Table** - Determines which quantification value is previewed in the **Create** table for SWE Segmentations.
 - **SWE Representation** – Sets which representation to display for SWE statistics: Shear Wave Velocity (m/s), Shear Modulus (kPa), or Young's Modulus (kPa).
 - **Default SWE Confidence Floor** – Default value for the confidence floor used to threshold SWE acquisitions upon load of SWE volumes.



- **Utilities and Display Options Settings** – Settings related to the image processing tools and display.
 - **Default Volume Projection** – Sets the default mathematical operation that will be used for volume projections.
 - **Show Mean Filter on Display (3D Static Only)** – When checked, automatically computes and activates the mean filter on loaded datasets.
 - **Default Mean Filter Radius** – Sets the default radius of the mean filter.
 - **Default AA Colormap** – Sets the default colormap for AA data.
 - **Default SWE Color Scalar Bar** – Sets the default SWE color scalar bar maximum and minimum values or (when checked) auto-scale maximum and minimum values for the selected SWE representation.
 - **Default Volume Renderer** – Sets the default option for 3D volume rendering.
 - **Default Rendering Technique** – Sets the default 3D rendering technique.
 - **Auto-scale Color Scalar Bar Range for SWE** – Sets whether the display range for SWE captures is automatically calculated or if the same set range is used for all SWE acquisitions.
 - **Display Volume Info on Screenshots** – When checked, volume information and metadata will be copied into screenshots and written into the output PNG file.
 - **Display Markup Statistics on Screenshots** – When checked, markup statistics will be copied into screenshots and written into the output PNG file.

6 Exporting Data

The following section describes different ways for exporting image data as well statistics produced in analysis.


6.1 Export Statistics

Exporting statistics creates spreadsheet files which include data from multiple markups. Each spreadsheet is for a different markup type (e.g. one for 3D segmentations, one for SWE segmentations). A spreadsheet can include markups from multiple volumes. The statistics for the spreadsheets can either be pulled from the scene, or from statistics saved to file. Statistics in the scene exist in Quantify tables for loaded volumes, and in caliper lengths. Statistics saved to file exist in markup folders, with the name “stats.csv”.

6.1.1 Export Statistics from Scene

Statistics can be exported from the scene in two main ways:

1) Clicking the **Export All Statistics** button at the bottom of the Quantify tab will export statistics for all loaded volumes that have markups or calipers, not only the primary and secondary volumes.



The screenshot shows the 'Quantify' tab in the SonoEQ software. It displays two sections: 'Surface statistics' and 'Image statistics (Primary)'. Both sections show data for 'Kidney_left' and 'Kidney_right' volumes. The 'Surface statistics' table includes columns for Voxels, Volume (mm^3), Volume (cc), and Surface Area (mm^2). The 'Image statistics (Primary)' table includes columns for Positivity (%), Mean (count), Median (count), and Standard Deviation (count). Below the primary statistics, there is a section for 'Image statistics (Secondary)' with similar columns. At the bottom of the Quantify tab, a green button labeled 'Export All Statistics' is highlighted with a red dashed box.

Voxels	Volume (mm ³)	Volume (cc)	Surface Area (mm ²)	
Kidney_left	1.65E+5	171.505	0.172	160.689
Kidney_right	2.42E+5	251.967	0.252	221.349

Positivity (%)	Mean (count)	Median (count)	Standard Deviation (count)	
Kidney_left	100.000	46.830	47	13.737
Kidney_right	100.000	52.636	52	13.460


Positivity (%)	Mean (count)	Median (count)	Standard Deviation (count)	
Kidney_left	100.000	52.173	52	20.174
Kidney_right	100.000	65.621	65	18.842

Export All Statistics

Clicking the Export All Statistics button will prompt a save file dialog where you can select where to save the data. Three comma-separated value (CSV) files will be saved:

- segmentation_stats_YYYY_MM_HH_MM_SS.csv - 3D segmentation stats.
- caliper_stats_YYYY_MM_HH_MM_SS.csv - caliper stats.
- cardiac_stats_YYYY_MM_HH_MM_SS.csv - M-Mode cardiac curves and 3D segmentation sequences.
- swe_stats_YYYY_MM_HH_MM_SS.csv – Shear Wave Elastography Statistics from SWE segmentations


2) In the File menu on the top toolbar, you can select to export all statistics from the **Export Statistics from Scene** Menu (same as clicking the Export All Statistics button) or individual types:



6.1.2 Export Statistics from Folder

Exporting statistics from saved files can be done in the following way:

- 1) Click File → **Export Statistics** from Folder



- 2) A data dialog will popup that behaves the same way as the Load Dialog; using **Select Folder** and **Select File(s)** will add statistic files (“stat.csv” files) to the dialog, which can be filtered.

Select Folder

Select File(s)

Filter...

		File Paths	Statistic Type
1	<input checked="" type="checkbox"/>	Test Study - TP1 - animal1 - [10:15:46] - M-Mode - Live - Markup_1	Cardiac
2	<input checked="" type="checkbox"/>	Test Study - TP1 - animal1 - [10:20:07] - SWE - 2D - Markup_2	SWE Segmentations
3	<input checked="" type="checkbox"/>	Test Study - TP1 - animal1 - [11:13:51] - AA - 3D - Markup_1	3D Segmentations
4	<input checked="" type="checkbox"/>	Test Study - TP1 - animal1 - [11:06:49] - B-Mode - 3D - Corrected: Medium - C_1	Caliper
5	<input checked="" type="checkbox"/>	Test Study - TP1 - animal1 - [10:23:44] - B-Mode - 3Dt - Phase Gated - Markup_1	Cardiac

Clear Table

Select Statistic Type:

All

All

3D Segmentations

2D Segmentations

Cardiac

Caliper

SWE Segmentations

File Count: 6

Export

Close

- 3) Click the drop-down list next to “Select Statistic Type:” to specify what type of statistics you wish to export. Only files visible in the table will be included in the export. You can also uncheck specific statistic files you want to exclude from the export. When ready click the **Export** button in the lower right of the dialog.

Copyright © 2023 Revvity, Inc.

6.2 Exporting Image Data

Image data can be exported to a new location inside SonoEQ. Clicking **Export Data** in the File menu on the top tool bar will bring up a window where specific data can be selected for export:

Data from any subfolder can be chosen and the file structure will be maintained when exporting so it can be loaded into SonoEQ. Checking the **Include Raw Files in Export** setting will include raw files used for reconstructing the final volumes. Checking the **Send to Compressed (zipped) folder** will export the data to a zip folder.




NOTE: Some of the .MHD files will be exported even if the **Include Raw Files in Export** option is not selected because they contain import information about the volumes.

6.3 Study Report Generator

SonoEQ can generate a PDF report of a study which tracks a markup statistic over time. To generate a report:

- 1) In the file menu, click **File → Create Study Report**.



- 2) A Study Report Generator window will pop up. In the Study Name Folder field click the (...) button and select a study folder. This folder's name will be the same as the **Study Name** field in acquisition.

 A screenshot of the 'Study Report Generator' dialog box. It contains several input fields: 'Study Name Folder' (empty, with a dropdown arrow and an ellipsis button), 'Create Report In' (set to 'C:\Users\VegaUser\Documents', with a dropdown arrow and an ellipsis button), 'Imaging Mode' (set to 'B-Mode'), 'Dimensionality' (set to '3D Static'), and 'Statistic' (set to 'Volume (mm^3)'). Below these fields is a section labeled 'Plot selected markups over time:' followed by a large empty rectangular box. At the bottom right are 'Generate' and 'Close' buttons. A red arrow points to the ellipsis button next to the 'Study Name Folder' field.

- 3) The box labeled "Plot selected segmentations over time:" will become filled with every markup name that exists for the selected Imaging Type and Dimensionality.

Study Name Folder: C:\Users\VegaUser\Documents\SonoEQ\Study Name

Create Report In: C:\Users\VegaUser\Documents

Imaging Mode: B-Mode

Dimensionality: 3D Static

Statistic: Volume (mm³)

Plot selected markups over time:

<input type="checkbox"/>	Kidney
<input type="checkbox"/>	Liver

Generate Close

- 4) Click the (...) button for the “Create Report In” text box and choose where you want the PDF file to be saved.

Create Report In: C:\Users\VegaUser\Documents


- 5) Select the desired **Imaging Type** from the drop-down menu. Then select the desired **Dimensionality** from the drop-down menu. Then select the desired **Statistic** to be in the report.
- 6) Now select which markups you want to be in the report by clicking the checkboxes in the “Plot selected segmentations over time:” table.

Plot selected markups over time:

<input type="checkbox"/>	Kidney
<input type="checkbox"/>	Liver

Generate Close


- 7) Click the **Generate** button and your PDF report will be generated in the selected **Create Report In** folder. Each page of the PDF will show an image of the markup, a plot of the selected statistic over time, and a list of details about the scan and markup on that page.



Study Name: Ortho Pancreas
Timepoint: 2020_11_04
Animal ID: 163
Acquired At Time: 2020_11_04-10_02_40
Imaging Mode: B-Mode
Dimensionality: 3D Static
Days Since First Acquisition: 7.1
Markup Name: Markup_1
Volume (mm³): 52.6


All Markup Statistics:

Voxels:	135994.0
Volume (mm ³):	52.56
Volume (cc):	0.05
Surface Area (mm ²):	71.61
Positivity (%):	100.0
Mean (count):	62.14
Median (count):	59.0
Standard Deviation (count):	24.19
Max (count):	213.0
Min (count):	2.0



7 Tools


The tools menu has several options which will be explained here.



7.1 Data Processing


The first 3 menu items include tools for reprocessing data, which work in a similar manner as the Load Dialog; a dialog is populated using **Select Folder** and **Select File(s)**, the files to re-process can be fine-tuned by filtering, and checking/unchecking entries in the table and clicking **Run** will reprocess the files selected on the dialog. These tools can be used to process volumes for the first time (if Post-Scan setting “Fuse Options” is set to “Fuse Later” and you choose not to process data upon closing SonoEQ) or reprocess data with updated algorithms.

Volume Reconstructor – reconstructs 3D acquisitions from raw files using calibration files saved on system.



- **Threshold for Breathing Motion** – remove frames corrupted by breathing. Can be set to “Off”, “Low”, “Medium”, “High”. The threshold for the algorithm defines how aggressively to remove artifacts with “High” being the most aggressive, which might result in removing data without breath artifacts.

Run Phase Gating – phase corrects cardiac 2D or 3D B-Mode sequences. Parameters can be changed individually for each entry on the table or the values of entry currently selected can be applied to all entries by clicking **Apply to all files**.



Select Folder		Select File(s)	
Filter...			
<input checked="" type="checkbox"/>	Volume		
1	<input type="checkbox"/> Test Study - TP2 - animal 1 - [10:28:15] - B-Mode - 2Dt	Processing Parameters Test Study - TP2 - animal 1 - [10:53:35] - B-Mode - 2Dt Cardiac Motion Freq Range (BPM) <input type="text" value="310.0,750.0"/> Auto Hodrick-Prescott Filter Smoothing <input checked="" type="checkbox"/> Hodrick-Prescott Filter Smoothing <input type="text" value="17"/> Axial Downsample Factor <input type="text" value="3"/> Axial Downsample Factor Phase <input type="text" value="8"/> Frame Rate Factor <input type="text" value="2"/> <input type="button" value="Apply to all files"/>	
2	<input checked="" type="checkbox"/> Test Study - TP2 - animal 1 - [10:28:43] - B-Mode - 3Dt		
3	<input checked="" type="checkbox"/> Test Study - TP2 - animal 1 - [10:53:35] - B-Mode - 2Dt		
<input type="button" value="Clear Table"/> File Count: 1		Total File Size: 13.0 MB <input type="button" value="Run"/> <input type="button" value="Cancel"/>	

- **Cardiac Motion Freq Range (BPM)** – Range of heart rates to be analyzed in beats per minute.
- **Auto Hodrick-Prescott Filter Smoothing** – Automatically calculate Hodrick-Prescott filter smoothing value.
- **Hodrick-Prescott Filter Smoothing** – Coefficient for separating breaths from cardiac motion. Higher values generally work for slower heart rates, or **Frame Rate Factor** less than 2. This parameter is used only if **Auto Hodrick-Prescott Filter Smoothing** is not checked.
- **Axial Downsample Factor** – Amount of spatial downsampling of original data for making single cardiac cycle sequence. A lower value may result on better image quality but will increase the processing time.
- **Axial Downsample Factor Phase** – Amount of spatial downsampling for separating breaths from cardiac motion. A lower value may result in better separation but will increase processing time.
- **Frame Rate Factor** – amount of temporal downsampling of original sequence. Lower value may result in better separation of breaths from cardiac motion but will increase processing time.

Run SWE Processing – runs SWE processing on 2D or 3D SWE raw data.

Select Folder

Select File(s)

Filter...

			Processing Parameters
	<input checked="" type="checkbox"/>	Volume	
1	<input checked="" type="checkbox"/>	Test Study - TP2 - animal1 - [08:47:13] - SWE - 2D	Test Study - TP2 - animal1 - [08:45:48] - SWE - 3D
2	<input type="checkbox"/>	Test Study - TP2 - animal1 - [08:49:15] - SWE - 2D	<div>Push Offset (mm) <input type="text" value="0.8"/></div> <div>Push from side <input type="text" value="Left"/></div>
3	<input checked="" type="checkbox"/>	Test Study - TP2 - animal1 - [08:45:48] - SWE - 3D	

Clear Table

File Count: 2

Reset to file metadata

Apply to all files

Total File Size: 120.0 MB

Run

Cancel

- **Push Offset (mm)** – Distance between push and tracking region.
- **Push From Side** – Determines on which side the push is in relation to the tracking region.

“Push Offset (mm)” and “Push from side” are enabled for volumes collected using a version of SonoEQ older than v1.12.0 and can be set the values of the current entry for all entries by clicking **Apply to all files**.


7.2 Volume Compressor

Clicking Volume Compressor opens a dialog that can be populated using **Select Folder** and **Select File(s)** buttons, like for the Load Dialog and data processing dialogs, for compressing and uncompressing raw image data to save disk space.


Select Folder		Select File(s)	
Filter...			
<input checked="" type="checkbox"/>	Volume		
1	<input checked="" type="checkbox"/>	Test Study - TP2 - animal 1 - [15:28:18] - AA - 3D	
2	<input checked="" type="checkbox"/>	Test Study - TP2 - animal 1 - [08:57:50] - B-Mode - 2D	
3	<input checked="" type="checkbox"/>	Test Study - TP2 - animal 1 - [08:57:50] - SWE - 2D	
Clear Table		Compress ▼	File Count: 3 Total File Size: 240.4 MB Run Cancel

8 Instrument

The **Instrument** menu has several options which will be explained here.



Vega without a heater



Vega with a heater

8.1 Servicing Dashboard

The **Servicing Dashboard** is used to service Vega systems and is available only to Revvity technicians.

8.2 Heat Imaging Bays

On Vega instruments configured with the internal heater module, an option will be present to **Heat Imaging Bays**. Clicking this menu option will turn the heater ON or OFF.



CAUTION:

Turning the heater OFF when an animal is sedated on the Vega can lead to hypothermia and harm the animal. Use caution when turning the heater OFF and ensure animal homeostasis is maintained via other methods.

8.3 Initialize Motion Stage

Selecting this menu option will start the motor system initialization homing procedure. It is only available when **Image Acquisition** connections have been successfully established and can be used if manual re-homing of the instrument is desired.

8.4 Firmware Updater

This feature allows updating or reverting the firmware of the wobbler engine, linear array engine, and heater controller.



NOTE:

If your "Observed Firmware" is the latest available, and you would like to revert back to a previous version of SonoEQ (e.g. pre-v1.13) you will need to revert back to the older firmware before SonoEQ can be used.

Wobbler Engine

Observed Firmware: ?

Linear Array Engine

Observed Firmware: ?

Heated Imaging Bays Controller

Configuration Bits: ?

Observed Firmware: ?

Firmware outdated

Wobbler Engine

Observed Firmware: ✓

Linear Array Engine

Observed Firmware: ✓

Heated Imaging Bays Controller

Configuration Bits: ✓

Observed Firmware: ✓

Firmware check OK
:03000000FFDAFD27
:000000001FF

Firmware current

- **Check** – Clicking this button will update the corresponding icon to report the status of the firmware. Hovering over this icon after checking will display a tooltip with information about the current firmware.
- **Observed Firmware** – Found by hovering over the status icon after clicking the “check” button. This is the version of the firmware currently installed on the system.
- **Latest Available Firmware** – Found by hovering over the status icon of the linear array engine after clicking the “check” button. This is the most up to date version of the firmware. If the “Latest Available” and “Current Flashed” do not match, updating may be required to image with your system.
- **Revert to pre-1.13.0 Version** – Clicking this button allows you to roll back the firmware version and install the compatible firmware for versions of SonoEQ prior to 1.13.0.
- **Update to latest available** – Clicking this button initiates the wizard that will guide you through the firmware upgrade process.

revvity	SonoEQ Version: 2.0.2		Revision: 22
	Title: SonoEQ User Guide		
	Proprietary & Confidential	Effective Date: 11/22/2023	Page 103 of 108

8.4.1 To revert firmware for use with pre-1.13 SonoEQ:

1. Open SonoEQ 1.13+ (only SonoEQ 1.13+ has the Firmware Updater)
2. Click Instrument → **Firmware Updater**
3. Click “Revert to pre-1.13.0 version”
4. When prompted, turn off the Vega for 30 seconds, then turn in back on
5. Re-open the Firmware Updater and check that the “Observed Firmware” has an earlier date than the “Latest Available Version”
6. Close SonoEQ 1.13+ and install the pre-1.13 SonoEQ version you would like to use.

9 Data Storage

Unless otherwise specified, all data captured with SonoEQ is stored under:

%UserProfile%\Documents\SonoEQ

We recommend regular backups of the data files to ensure the best user experience when running SonoEQ.

If available, it is also recommended to backup data to a network-attached storage (NAS) server or cloud-based storage.

File Format Reference

This list defines the various file formats that are produced by SonoEQ.

Extension	Description
.mha	Binary file with reconstructed 3D image data
.raw	Binary file with raw, non-reconstructed 3D image data
.mhd	Metaheader data file associated with a .raw file, necessary for reconstruction
.seq.*	Metaheader data file of sequence data (i.e. CINE, 4D). The * represents either .mhd or .nrrd (e.g. ".seq.mhd")
.seg.nrrd	Nearly-raw raster data (NRRD) file that contains 3D volume data used for a single 3D segmentation
.csv	Comma-separated value file for exported statistics
.txt	Text file with saved study notes
.png, .tif	Image files for thumbnail projections of reconstructed 3D image data, and screenshots of camera view and ultrasound stream



CAUTION:

Be careful when moving or renaming the SonoEQ folder structure using Windows Explorer. SonoEQ relies on the folder structure to load associated files (i.e. a segmentation that was done on a specific image volume), and changes to this folder structure may result in unexpected outcomes.

revvity	SonoEQ Version: 2.0.2		Revision: 22
	Title: SonoEQ User Guide		
	Proprietary & Confidential	Effective Date: 11/22/2023	Page 105 of 108

10 Mouse and Keyboard Shortcuts

Below is basic information about how to use the three-button mouse on Windows to perform operations in SonoEQ. If a specific *View* is indicated in the table below, the corresponding view must be clicked on first before the keyboard shortcut will function.

Mouse Operations

Action	View	Mouse sequence
Rotate	3D-only	Left-Click & drag OR
	Slice-only	Ctrl + Left-Click & drag (Rotate in-plane)
Zoom	All views	Right-Click & vertical drag
Pan	All views	Middle-Click & drag OR
		Shift + Left-Click & drag
Brightness/ Contrast	All views	Press Alt + “c” to enter adjustment mode then Left-Click & drag <ul style="list-style-type: none"> Adjust the brightness by moving the mouse vertically Adjust the contrast by moving the mouse horizontally
Cross Reference	Slice-only	Hold “Shift” while moving the mouse in any Slice Viewer will cause other Slice Viewers to scroll to the same position
Maximize View	All views	Double Left-Click to maximize slice viewer. Double Left-Click again to minimize.
Copy Image	All views	Right-Click on slice viewer and Left-Click on Copy Image to copy image to palette
Expand ROI	Slice-only	Press Alt + Left-Click & drag ROI points to expand/contract ROI

revvity	SonoEQ Version: 2.0.2		Revision: 22
	Title: SonoEQ User Guide		
	Proprietary & Confidential	Effective Date: 11/22/2023	Page 106 of 108

Keyboard Operations

Key	View	Effect
r	All views	Fits the camera's field of view to show all that is visible in the viewer
v	Slice-only	Toggles the slice plane visibility in the 3D view
h	Slice-only	Cycle through cardiac curves
c	All views	Add Caliper
Insert	All views	Create new markup
Delete	All views	Delete selected markup


Creating and Editing Segmentations

Key	View	Effect
a	All views	Apply segmentation
e	All views	Edit segmentation
g	All views	Toggle active segmentation visibility
z	All views	Delete last control point placed
Space Bar	All views	Toggle Placement mode
Esc	All views	Cancel changes to segmentation

revvity	SonoEQ Version: 2.0.2		Revision: 22
	Title: SonoEQ User Guide		
	Proprietary & Confidential	Effective Date: 11/22/2023	Page 107 of 108

11 Troubleshooting

In case of software crashes or unexpected behavior, you can access SonoEQ's log files to share with the Customer Support team (<http://www.revvity.com/contact-us>) by using the **Report a bug** feature, which can be accessed through **Help** on the top toolbar.



The most recent logs can be viewed in the “Recent log files” section of the “Report bugs and request enhancements” window. The contents of a selected log are displayed in the “Log file content” section. The **Copy log to clipboard** will allow you to copy and paste the log file when emailing Customer Support. The **Open log file location** will open a Windows explorer window at the location where the log files are saved so you can attach logs when emailing Customer

revvity	SonoEQ Version: 2.0.2		Revision: 22
	Title: SonoEQ User Guide		
	Proprietary & Confidential	Effective Date: 11/22/2023	Page 108 of 108

Support. The **Open log in file editor** opens the currently selected log in a text editor such as Notepad so you can save it in a location of their choice.